
PHASITRON®

Breathing circuit for VDR®-4

EN FR IT NL

PERCUSSIONAIRE®
CORPORATION

130 McGhee Rd., Suite 109, Sandpoint, Idaho 83864
MDSS GmbH, Schiffgraben 41, 30175 Hannover, Germany
percussionaire.com +1.208.263.2549
©2020 Percussionaire® Corporation P20067 Rev A
All rights reserved.

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Intended Use
The Phasitron® breathing circuit kit is intended to be used for continuous ventilation intended for controlled ventilation of patients.

Patient Population
Neonates, pediatric, and adults.

Indications for Use
The Phasitron® breathing circuit kit is intended to be used with Percussionaire® control device model VDR®-4. Indicated for patients diagnosed with either respiratory failure or respiratory insufficiency requiring partial or full ventilator support in hospitals.

Absolute Contraindications:
• Untreated tension pneumothorax
• Untrained or unskilled operator

Relative Contraindications:
• History of pneumothorax • Pulmonary air leak
• Recent pneumectomy (without functioning chest tube)
• Pulmonary hemorrhage • Lack of patient cooperation
• Myocardial Infarction • Vomiting

Possible Adverse Reactions
• Hyperventilation • Gastric distension
• Decreased cardiac output • Increased intracranial pressure
• Increased air trapping • Contamination
• Hyperoxygenation • Pneumothorax
• Pulmonary air leak • Pulmonary hemorrhage

Document Symbols

	WARNING		Type BF Applied Part
	CAUTION		Single Patient Use
	Read the manual before use		Prescription Only
	CE marking		Catalog Number
	Manufacturer		Lot Number
	Manufacturer Date		European Representative
	Non-Sterile		Not Made with Natural Rubber Latex
	Does Not Contain the Phthalate Plasticizers DEHP, DBP, DBP or BBP		Disposal


Read and understand the instructions completely before setting up and using your Phasitron® breathing circuit kit.

Precautions


- CAUTION** icon indicates a risk of equipment damage.
- WARNING** icon indicates a risk of injury to patient or operator.
- CAUTION: Do Not autoclave or expose to processing temperatures above 50°C/122°F**
- WARNING: Do not occlude the exhalation or entrainment valve gate ports.**
- Type BF**
- Single patient use**
- All persons providing HFV must be trained in the use of Percussionaire® devices.
- Therapists should evaluate how their patients tolerate HFV. Auscultation and observation of the mechanical vibrations of the chest and abdomen are primary indicators of effective ventilation delivery.
- Only use genuine Percussionaire® parts and accessories.
- Never use external fixation to keep the Phasitron®, mask, or mouthpieces on the patient's face when using in noninvasive modes.

Phasitron® Setup Connecting to the VDR®-4


- Place heater chamber onto heater. Insert nebulizer into mounting hook, normally installed near heater.
- Connect red, clear, green, and yellow tubing quick-connect fittings to the VDR®-4 front panel (Fig. 1). (May come pre-assembled)


- Connect blue nebulizer connection, without the O-ring, to the heater chamber using unheated fresh gas blue corrugated tubing (not supplied with kit) (see Fig. 2, A and also see Fig. 3, C). Connect the blue heated wire tubing to heater (not supplied with kit) (Fig. 2, B).

- Insert 4-way manifold 22mm OD port into the O-ring side of the blue nebulizer placed upright (Fig. 3, D). (May come pre-assembled)

- Install 3L anesthesia bag 22mm ID onto bottom of manifold (Fig. 3, E).
- Replace 3L anesthesia bag with 0.5L anesthesia bag for neonatal to pediatric patients, if required.
- Connect white expiratory corrugated heated wire limb (not supplied with kit) to the side inlet port of manifold (Fig. 3, F).
- Connect green tubing quick-connect fitting to the top of blue nebulizer tee green fitting (Fig. 3, G).
- Connect yellow tubing quick-connect fitting to bottom of nebulizer bowl (Fig. 3, H).
- Remove wye from heated wire circuit (not included in kit) and attach white corrugated tubing to the clear expiratory port valve on the Phasitron® (Fig. 4, M).

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