
# PHASITRON®

## Breathing circuit A50605-D

EN DE FR IT NL

**PERCUSSIONAIRE®**  
CORPORATION

130 McGhee Road, Suite 109, Sandpoint, Idaho 83864  
percussionaire.com +1.208.263.2549

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### Intended Use

The Phasitron® breathing circuit kit is intended to be used for continuous ventilation for the controlled ventilation of patients.

### Patient Population

Neonates, pediatric, and adults.

### Indications for Use

The Phasitron® breathing circuit kit is intended to be used with a Percussionaire® HFPV system for either hospital or prehospital use where emergency care is being provided, including intrahospital or external hospital transport.

### Absolute Contraindications

- Untreated tension pneumothorax
- Untrained or unskilled operator

### Relative Contraindications

- History of pneumothorax
- Recent pneumonectomy
- Pulmonary hemorrhage
- Myocardial infarction
- Pulmonary air leak (without functioning chest tube)
- Lack of patient cooperation
- Vomiting

### Possible Adverse Reactions

- Hyperventilation
- Decreased cardiac output
- Increased air trapping
- Hyperoxygenation
- Pulmonary air leak
- Gastric distension
- Increased intracranial pressure
- Contamination
- Pneumothorax
- Pulmonary hemorrhage

### Document Symbols

⚠️ <b>WARNING</b>	📄 <b>Type BF Applied Part</b>
⚠️ <b>CAUTION</b>	👤 <b>Single Patient Use</b>
📖 <b>Read the manual before use</b>	📖 <b>Prescription Only</b>
📑 <b>CE marking</b>	📑 <b>Catalog Number</b>
🏢 <b>Manufacturer</b>	🏢 <b>Lot Number</b>
📅 <b>Manufacturer Date</b>	📅 <b>European Representative</b>
🧼 <b>Non-Sterile</b>	🚫 <b>Not Made with Natural Rubber Latex</b>
🗑️ <b>Does Not Contain the Phthalate Plasticizers DEHP, DBP, DIBP</b>	♻️ <b>Disposal</b>



📖 **Read and understand the instructions completely before setting up and using your Phasitron® breathing circuit kit.**

### Precautions

- ⚠️ WARNING: Indicates the possibility of injury to the user or operator.**
- ⚠️ CAUTION: Indicates the possibility of damage to the device. US Federal law restricts this device to sale by or on the order of a physician.**
- ⚠️ CAUTION: Do not expose to processing temperatures above 50 °C/122 °F.**
- ⚠️ WARNING: Single Patient Use, not intended for reprocessing.**
- ⚠️ WARNING: Do not occlude the exhalation or entrainment ports.**
- All persons providing mechanical ventilation must be trained in the use of Percussionaire® High-Frequency devices.
- Therapists should evaluate how their patients tolerate the treatment. Auscultation and observation of the mechanical vibrations of the chest and abdomen are primary indicators of effective ventilation.
- Only use genuine Percussionaire® parts and accessories.
- Never attempt to modify unit.

### Phasitron® Setup

Connect breathing circuit to the High-Frequency Percussive Transport ventilator.

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All rights reserved. P20122 Rev A