# In-Line Valve

## Instructions for Use

- EN
- JA
- RU

**Percussionaire®**  
130 McGhee Rd., Suite 109,  
Sandpoint, ID 83864 USA


For contact details of your local Sentec partner, visit [sentec.com/contact/](http://sentec.com/contact/)

© 2025 Percussionaire® All rights reserved.  
P20020 Rev H | 2025-09

# sentec.

The In-Line Valve includes a cap to close the Phasitron® Port.

### Treatment Frequency Using IPV® In-Line Therapy

IPV® in-line use with a ventilator is based on patient need. 2 times per day or up to every 4 hours is recommended. Always use institutional/hospital protocol when possible.

- ⚠ **WARNING:** Follow your institutional protocols for infection control.
- ⚠ **WARNING:** Follow institutional protocols before disconnecting ventilator inspiratory limb prior to installation of IPV® In-Line Valve.

### Initial Placement of IPV® In-Line Valve

- ⚠ **WARNING:** The In-Line Valve should not be used in neonates\* or infants\* who are at high risk of cardiopulmonary or neurological compromise due to unrecognized hypoventilation. \*Indications for the use of IPV therapy in patients under two years of age vary by geography and by device.
- ⚠ **CAUTION:** Ensure pressure relief adjustment valve is closed.
- ⚠ **CAUTION:** Ensure IPV® In-Line Valve is inserted into inspiratory side of ventilator circuit.

1. Install the IPV In-Line Valve into the inspiratory limb of the ventilator circuit in either of the following configurations:

A. As close to the patient wye as allowable

B. On the dry side of the humidifier

2. Insert the provided cap into the Phasitron® port of the In-Line Valve, and verify the pressure relief port is fully closed.
- NOTE: The IPV In-Line Valve may create a leak in the ventilator circuit, even if the pressure relief knob is fully closed, which may reduce delivered (exhaled) tidal volumes.
3. Allow the ventilator to complete two or more cycles and adjust settings, if necessary, to compensate for any potential leaks and to ensure exhaled tidal volumes are adequate.

⚠ **WARNING:** Insufficient exhaled tidal volumes may lead to hypoventilation.

NOTE: If adequate exhaled tidal volumes cannot be achieved with the In-Line Valve in place, consider removing the In-Line Valve from the circuit when IPV therapy is not being delivered or choose a different treatment configuration (e.g., direct to artificial airway).

### Adding Phasitron 5® to In-Line Valve


⚠ **WARNING:** NEVER perform IPV treatment without liquid in Phasitron® nebulizer during treatment.

1. Fill the Phasitron® nebulizer with 15 to 20 cc of saline solution or other liquids.

NOTE: Aerosol consumption: approximately 0.75 - 1 ml per minute.

2. Remove the cap from the In-Line Valve's Phasitron® port and insert the Phasitron®.

Installation of IPV® In-Line Valve with Phasitron® 5

⚠ **WARNING:** Expiratory port on Phasitron® 5 is to remain unobstructed throughout therapy when used in line or direct to airway.

### Delivering IPV® Therapy

- ⚠ **WARNING:** If chest percussion is inadequate, increase inspiratory flow or raise operational pressure and scan frequency rate to mobilize secretions.
- ⚠ **WARNING:** Source pressure, percussive frequency, and inspiratory flow can be adjusted to increase and decrease amount of "chest wiggle".
- ⚠ **WARNING:** Suctioning should be performed as needed. Note the current ventilator alarm and mode settings.
- ⚠ **WARNING:** Expiratory port on Phasitron® 5 is to remain unobstructed throughout therapy when used in line or direct to airway.
- ⚠ **CAUTION:** Note the current ventilator alarm and mode settings. Pressure-control mode is recommended or follow institutional protocol.
- ⚠ **CAUTION:** Mean Airway Pressures will increase slightly with the administration of IPV® in-line therapy with the ventilator. The respiratory care practitioner needs to be aware of this effect and monitor the patient closely for any adverse side effects.
- ⚠ **CAUTION:** When using the ventilator in pressure-control mode, the In-Line Valve may remain closed. When using the ventilator in volume-control mode, the In-Line Valve may be opened to create a leak.
- ⚠ **CAUTION:** May take multiple treatments to identify optimal therapeutic effect for each patient.
- ⚠ **CAUTION:** Reset occasional high-pressure alarms as they occur. When applying IPV® in line, adjust the pressure relief adjustment

(continued on other side)

## EN IPV® In-Line Valve, 22 mm ID x 22 mm OD Single Patient, Instructions for Use

### Intended Use

The IPV® In-Line Valve is intended to be used to provide IPV® (Intrapulmonary Percussive Ventilation) therapy to patients while assisted by Conventional Mechanical Ventilation (CMV) using pressure-control, volume-control, SIMV-PC etc., when direct connection of IPV® is not indicated.

For use only with Percussionaire® IPV® Control Devices. Compatible with all Percussionaire® single-patient Phasitrons.

### Patient Population

Neonates to adults, for whom IPV® therapy has been prescribed.

### Absolute Contraindications:

- Untreated tension pneumothorax
- Untrained operator

### Relative Contraindications:

- History of pneumothorax
- Myocardial Infarction
- Vomiting
- Recent pneumonectomy
- Pulmonary hemorrhage
- Pulmonary air leak (Without functioning chest tube)

### Possible Adverse Reactions:

- Decreased cardiac output
- Pneumothorax
- Hyper-oxygenation
- Pulmonary hemorrhage
- Gastric distension
- Increased intracranial pressure
- Increased air trapping
- Pulmonary air leak
- Hyperventilation

### Symbols

⚠ <b>WARNING</b>	⚠ <b>Type BF Applied Part</b>
⚠ <b>CAUTION</b>	⚠ <b>Single Patient Multiple Use</b>
📖 <b>Read the manual before use</b>	📖 <b>Prescription Only</b>
🏢 <b>Manufacturer</b>	📄 <b>Catalog Number</b>
📅 <b>Manufacture Date</b>	📄 <b>Lot Number</b>
	🕒 <b>Use-by Date</b>

📖 Read these instructions before setting up and using your IPV® In-Line Valve kit.

### Precautions


- ⚠ **WARNING:** Indicates the possibility of injury to the user or operator.
- ⚠ **CAUTION:** Indicates the possibility of damage to the device.

US Federal law restricts this device to sale by or on the order of a physician.

- All persons providing IPV® treatment must be trained in the use of Percussionaire®, Inc. devices.
- Therapist should evaluate patient's response and tolerance to therapy. Adjust inspiratory flow or operating pressures up or down as appropriate. Auscultation and observation of the mechanical vibrations of the chest and abdomen are primary indicators of effective treatment. Utilize ETCO<sub>2</sub> before and after treatments and SaO<sub>2</sub> as secondary indicators.
- Use only genuine Percussionaire®, Inc. parts and accessories.

### P5-TEE In-Line Valve

Pressure relief adjustment valve

valve to achieve desired institutional/hospital protocol Peak Inspiratory Pressure (PIP). High-pressure alarms should not occur on a regular basis if the Pressure relief adjustment valve is set correctly.

Consult the IFU of your IPV® device for instructions on how to deliver IPV® Therapy.

NOTE: Patients who are performing T-tube trials or CPAP sprinting may be taken off of the ventilator for the IPV® treatment utilizing a flex adapter. Decreasing cuff pressure still applies to this patient population.

NOTE: Following your institutional protocols for cuffed endotracheal tubed patient, the cuff pressure may be lowered.

NOTE: Lowering of the cuff pressure facilitates secretion removal into the oral cavity where they may be suctioned. This also helps in the prevention of tube obstruction in the event copious secretions are mobilized.

#### At Completion of Therapy

1. If cuff was deflated during treatment, reset cuff pressure.
2. Close pressure relief adjustment valve and turn off IPV device in quick succession.  
NOTE: The IPV In-Line Valve may create a leak in the ventilator circuit, even if the pressure relief knob is fully closed, which may reduce delivered (exhaled) tidal volumes.
3. Restore ventilator settings noted before treatment and reassess the patient.
4. Allow the ventilator to complete two or more cycles and adjust settings, if necessary, to ensure exhaled tidal volumes are adequate.

#### ⚠ WARNING: Insufficient exhaled tidal volumes may lead to hypoventilation.

NOTE: If adequate exhaled tidal volumes cannot be achieved with the In-Line Valve in place, consider removing the In-Line Valve from the circuit when IPV therapy is not being delivered or choose a different treatment configuration (e.g., direct to artificial airway).

5. Disconnect the Phasitron® from the In-Line Valve and re-insert the cap into the Phasitron® port of the In-Line Valve.

#### Specifications

Storage and Transport Range	Temp., -40°C to 5°C (-40°F to 41°F) without relative humidity control Temp., 5°C to 35°C (41°F to 95°F) with relative humidity up to 90% non-condensing
Service Life	30 days from date of opening Remove and dispose with ventilator circuit Single ventilator circuit use only
Shelf Life	2 years
Disposal	Recycle according to local laws.

#### Limited Warranty

The manufacturer warrants to the initial purchaser that each new In-Line Valve will be free from defects in workmanship and materials. The manufacturer's sole obligation under this warranty is to, at its own choice, repair or replace any component – for which the manufacturer acknowledges the warranty coverage – with a replacement component.

#### Warranty Exclusions and System Performance

Sentec can neither guarantee nor verify product performance characteristics, nor accept warranty claims or product liability claims, if the recommended procedures are not carried out, if the product has been subject to misuse, neglect, or accident, if the product has been damaged by extraneous causes, or if accessories other than those recommended by Sentec are used.

Any serious incident that has occurred in relation to the IPV® system has to be reported to Sentec (regulatory.percussionaire@sentec.com) and/ or to the competent authority of the country where the incident occurred. If unsure whether an incident is a reportable event, contact Sentec.

Percussionaire®, IPV® and Phasitron® are registered trademarks of Percussionaire®, Inc.

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