



PERCUSSIONAIRE®

# TXP<sup>5</sup> System

## User Manual



**High-Frequency Percussive Ventilator (HFPV)**



©2020 Percussionaire® Corporation

ALL RIGHTS RESERVED

1st Edition

First Printing May 2020

Percussionaire® is a registered trademark of Percussionaire Corporation

This work is the sole property of Percussionaire® Corporation. The information in this manual is confidential and may not be disclosed to third parties without the prior written consent of Percussionaire®. No part of this document may be copied, reproduced, transmitted, or stored in any electronic information system without the prior written permission of Percussionaire® Corporation.

One or more patents may cover the devices and products contained in this manual.

This manual was originally released and supplied in English.

For a list of available translations, contact [customerservice@percussionaire.com](mailto:customerservice@percussionaire.com)

All ventilators should be operated and serviced only by trained professionals.

Percussionaire® Corporation's sole responsibility with respect to its ventilators, accessories, components, and software, and their use, are as stated in the warranty provided in the manuals. The information set forth herein is believed to be accurate; it is not a substitute for the exercise of professional judgment.

## Table of Contents

<b>Chapter 1: Introduction</b>	1
High-Frequency Percussive Ventilation	1
Phasitron® Breathing Circuit A50606-TXP	1
Effects of HFPV	2
Document Symbols	3
Warnings and Cautions	4
⚠ Warnings	4
⚠ Cautions	4
<b>Chapter 2: Intended Use</b>	5
Indications for Use	5
Patient Population	5
Absolute Contraindications	5
Relative Contraindications	5
Possible Adverse Reactions	5
Physiological Benefits of HFPV	5
Clinical Limitations/Restrictions	6
<b>Chapter 3: System Description</b>	7
Frequency Control Knob	7
Amplitude Control Knob	7
Control Knob	8
Tubing Connectors	8
Digital Display	8
Power-On Self-Test (POST) Mode	8
Treatment Screen	9
Main Menu Options	9
Power Off	10
System Information	10
Contrast	10
Exit Menu	11
Volume	11
<b>Chapter 4: Alarms</b>	12
Alarm Settings	12
Alarm Conditions	14
High Priority	14
Battery Empty	14
Loss of Power Alarm	14
Medium Priority	15
Battery Low	15
Patient Pressure Below Set Limit	15
Alarm Silence	15
Battery Discharge	16
<b>Chapter 5: Setup</b>	17
Controller and Stand	17
DISS Gas Connection	17
Power Supply	17
Mounting to Pole/Post	18
Adjusting Clamp Range and Tension	18
Mounting to Rail	18

Phasitron® A50606-TXP Assembly .....	19
Connect Tubing Harness .....	19
Configurations .....	20
Intubated Patient with Phasitron® .....	20
Standard Circuit .....	20
Optional PEEP Valve Kit (PRT-A70143) .....	20
Optional Secondary Filters .....	20
<b>Chapter 6: Pre-Use Check</b> .....	21
<b>Chapter 7: Ventilation Settings</b> .....	22
Standard Circuit Initial Settings .....	22
Adjustment Options: Standard Circuit .....	22
Increase pO <sub>2</sub> and decrease CO <sub>2</sub> .....	22
Increase O <sub>2</sub> if CO <sub>2</sub> is OK .....	22
Decrease CO <sub>2</sub> if O <sub>2</sub> is OK .....	22
Adjustment Options: Standard Circuit with	
Optional PEEP Valve and Inspiratory Valve .....	23
Increase pO <sub>2</sub> and Decrease CO <sub>2</sub> .....	23
Increase O <sub>2</sub> if CO <sub>2</sub> is OK .....	23
Decrease CO <sub>2</sub> if O <sub>2</sub> is OK .....	23
Increase CO <sub>2</sub> if O <sub>2</sub> is OK .....	23
Patient Monitoring .....	24
Assessment .....	24
Ventilation .....	24
Increase Ventilation .....	24
Oxygenation .....	25
Increase Oxygenation .....	25
Decrease Ventilation: Increase pCO <sub>2</sub> .....	25
<b>Chapter 8: Cleaning and Maintenance</b> .....	26
Cleaning .....	26
TXP® 5 Controller .....	26
Digital Display Screen .....	26
Phasitron® A50606-TXP Breathing Circuit Kit .....	27
Cleaning and Disinfecting Solutions .....	27
Maintenance .....	27
TXP® 5 Controller .....	27
Lithium-ion Battery Care .....	27
Battery Replacement .....	28
<b>Chapter 9: Troubleshooting</b> .....	29
<b>Chapter 10: Technical Specifications</b> .....	31
TXP® 5 Ventilator .....	31
Measured Performance .....	31
FiO <sub>2</sub> Measurements .....	32
Digital Display Specifications .....	33
Phasitron® A50606-TXP .....	33
Battery .....	34
<b>Chapter 11: Service and Repair</b> .....	35
Disposal of Equipment .....	35
Battery Disposal .....	35
<b>Chapter 12: Limited Warranty</b> .....	35

## Chapter 1: Introduction

This chapter provides an overview of the TXP® 5 emergency high-frequency ventilator.




### High-Frequency Percussive Ventilation

High-Frequency Percussive Ventilation (HFPV) is a flow-regulated and time-cycled hybrid form of high-frequency ventilation which can be lifesaving for patients with fragile, small, congested, or stiff lungs weakened by barotrauma or infection. HFPV is a solution that requires neither the pressure of mechanical ventilation nor the breathing action of the patient. This HFPV system supports both diffusive and convective flow by stacking pulses in cumulative subtidal volumes, allowing for air exchange, airway clearance, and lung recruitment for most patients.


### Phasitron® Breathing Circuit A50606-TXP

HFPV provides respiratory assistance as the patient breathes through the Phasitron®. The patented Phasitron® uses a unique venturi mechanism to protect the lung from over-pressure. By automatically adjusting to the resistance of the lung, the Phasitron® precisely, and safely delivers the optimal amount and pressure of air required by the alveolar space. When lung resistance is low, as in a compliant lung, all the pulsed air from the TXP® 5 enters the mouth of the venturi. Each air pulse draws up to four times as much additional air into the venturi tube. This low-pressure entrained air automatically fills the available space in the lung. The Phasitron® continuously and instantaneously adjusts to keep a gentle and safe air pressure even in a compromised lung.

## Phasitron® Breathing Circuit A50606-TXP



















## Effects of HFPV



The effects of HFPV occur with or without the cooperation of the patient. HFPV provides a sub-tidal gas exchange within the respiratory bronchioles with associated alveolar recruitment maintaining a minimal mean intrathoracic expiratory pressure increase for peripheral lung stabilization. This allows for mechanical ventilation to provide for peripheral lung recruitment while minimizing the potential for induced barotrauma.

## Document Symbols

 <b>WARNING</b>	 <b>Type BF Applied Part</b>
 <b>CAUTION</b>	 <b>Single Patient Use</b>
 <b>Read the manual before use</b>	 <b>Prescription Only</b>
 <b>CE marking</b>	 <b>Catalog Number</b>
 <b>Manufacturer</b>	 <b>Lot Number</b>
 <b>Manufacture Date</b>	 <b>European Representative</b>
 <b>Non-Sterile</b>	 <b>Not Made with Natural Rubber Latex</b>
 <b>Does Not Contain the Phthalate Plasticizers DEHP, DIBP, DBP, or BBP</b>	 <b>Disposal</b>

## Warnings and Cautions

A **WARNING** icon indicates a risk of injury to a patient or operator.

A **CAUTION** icon indicates a risk of equipment damage.

### Warnings

<b>Clinical support</b>	If the HFPV device is used on a patient with an indwelling airway (i.e., endotracheal or tracheostomy tube), a clinician must be available at all times. These devices enhance secretion clearance. Patients must be assessed for a reduced functional residual capacity (FRC).
<b>Patient Monitoring</b>	Complete a clinical assessment before placing a patient on the TXP® 5 ventilator.
<b>Personnel Qualifications</b>	The operator of the TXP® 5 is responsible for reading and understanding the manual before use. TXP® 5 is a medical device designed for hospital use by trained clinicians under the supervision of a physician. Only change the prescription and device settings on the order of the supervising physician.
<b>Pre-Use Check</b>	Always conduct a pre-use check before starting HFPV with the TXP® 5.
<b>Suctioning</b>	Perform suctioning as necessary; pulmonary alveoli cannot be ventilated when their transmitting airways are obstructed.

### Cautions

<b>Cleaning</b>	Do not use any steam cleaning methods to clean the device or Phasitron® A50606-TXP breathing circuit. Always follow hospital/institutional protocols for cleaning and disinfection.
<b>Maintenance</b>	The TXP® 5 must not be opened by anyone other than Percussionaire®-authorized service personnel. Maintain and service the TXP® 5 device according to the recommendations provided in this manual. Only use Percussionaire® accessories explicitly designed for use with the TXP® 5 device.
<b>Malfunctions</b>	If the TXP® 5 malfunctions, it should not be utilized any further. Report any malfunctions immediately.
<b>Read and Understand Manual</b>	All persons providing HFPV therapy must be trained in the use of the TXP® 5.
<b>Safety</b>	Do not place objects on top of the TXP® 5. Do not cover the device during use. Do not lean on device.



## Chapter 2: Intended Use

### Indications for Use

The TXP® 5 is designed for emergency ventilation of adult, pediatric, or neonatal patients. The device is intended for hospital or pre-hospital use, emergency care, intra-hospital, and external hospital transport. The TXP® 5 is intended for use under the supervision of a licensed physician, by both properly trained clinicians and personnel with limited training. The TXP® 5 requires a 50 psi source capable of maintaining 15 LPM. A medical source of compressed oxygen or air is preferred, but other emergency sources may be used.

### Patient Population

TXP® 5 ventilator is intended for use on neonatal, pediatric, and adult patient populations.

### Absolute Contraindications

• Untreated tension pneumothorax	• Untrained or unskilled operator
----------------------------------	-----------------------------------

### Relative Contraindications

• History of pneumothorax	• Lack of patient cooperation
• Recent pneumonectomy	• Vomiting
• Pulmonary hemorrhage	• Pulmonary air leak (without functioning chest tube)
• Myocardial infarction	

### Possible Adverse Reactions


• Decreased cardiac output	• Increased intracranial pressure
• Pneumothorax	• Increased air trapping
• Hyper-oxygenation	• Pulmonary air leak
• Pulmonary hemorrhage	• Hyperventilation
• Gastric distension	


### Physiological Benefits of HFPV

• Recruitment of atelectatic lung	• Mechanical ventilation
• Improved FRC	• May improve breathing pattern
• Decreased work of breathing	• Increased secretion mobilization

## Clinical Limitations/Restrictions

Use of the TXP® 5 is limited to qualified individuals who have received training.

 **WARNING:** Suctioning should be performed as necessary; pulmonary alveoli cannot be ventilated when airways are obstructed.

 **WARNING:** When using an HME (Heat-Moisture Exchanger), connect between the Phasitron® breathing circuit and the patient. Follow all setup instructions provided by the HME manufacturer.


**NOTE:** *Any HME attached to the Phasitron® breathing circuit should comply with ISO 9360-1 or ISO 9360-2.*

## Chapter 3: System Description


The TXP® 5 ventilator provides high-frequency percussive ventilation (HFPV) as the patient breathes through the Phasitron® breathing circuit.



### Frequency Control Knob



### Amplitude Control Knob



## Control Knob



Use the control knob to select menus and silence alarms.


## Tubing Connectors



Connect the Phasitron® breathing circuit easily by inserting the tubing connectors into the bulkhead fittings.

**CAUTION:** Ensure the connection is straight to prevent crimping the O-ring and causing a leak.

## Digital Display



The treatment screen displays Amplitude, Mean Airway Pressure, Frequency, alarm high and low setpoints, battery condition, charging status, and power plug indicator.

The display has three menu options:

- Main menu
- Volume settings
- Alarm Settings

## Power-On Self-Test (POST) Mode


Push the control knob for more than .5 of a second to power on the display. The screen displays the POST screen while it evaluates system processes before proceeding to the treatment screen. The POST screen is displayed for 3 to 5 seconds.

**NOTE:** If a technical issue arises during POST, the System Failure screen will appear.

## Treatment Screen


Once POST is complete, the treatment screen is displayed.

The treatment screen displays several menu options. Turn the knob until the desired menu option highlights. Push the knob to select.




**NOTE:** Displayed values may read zero if the TXP® 5 is not on or pulsing.

## Main Menu Options




## Power Off



Highlight **POWER OFF** by turning the control knob. Push the control knob to select.

## System Information




Turn the control knob to highlight **SYSTEM INFO**. Push the control knob to select.




A popup will appear displaying system information. Turn the control knob to highlight **EXIT** to return to the treatment screen. Push the control knob to make selection.

## Contrast




Turn the control knob to highlight **CONTRAST**. Push control knob to select. Rotate control knob clockwise to increase (darken) or counterclockwise to lower (lighten).



Push control knob to **SET**.

## Exit Menu



Turn the control knob to highlight **EXIT**. Push the control knob to select and return to the treatment screen.

## Volume



Turn the control knob to highlight the volume icon on the treatment screen. Push the control knob to select.



A popup will appear displaying volume options. Turn the control knob to highlight the desired volume. Push control knob to select:

**LOW** – volume set to 25% of speaker output  
**MED** – volume set to 50% of speaker output  
**HIGH** – volume set to 100% of speaker output  
A confirming tone will sound when the volume selection is made.



Turn the control knob to highlight **EXIT**. Push to select. This will return the user to the treatment screen.


# Chapter 4: Alarms

## Alarm Settings


**NOTE:** Alarms are controlled by the Mean Airway Pressure settings.




Turn the control knob to highlight the alarm icon. Push the control knob to select.



Alarm default set to:  
High 45  
Low 3  
An alarm popup will appear.






Turn the control knob to highlight the value to be changed. Push to select.



Turn the control knob to raise or lower the alarm value. Push to select.



	<p>When satisfied with the new alarm settings, turn the control knob to highlight the <b>SAVE</b> button. Push to select. This will return the user to the treatment screen, or select <b>QUIT</b> to return to the treatment screen without updating or changing alarm settings.</p>
	<p>Highlight and select <b>EXIT</b> to return to the treatment screen without updating or changing alarm settings.</p>
	<p>The alarm parameters display on the treatment screen.</p>


**NOTE:** The process is the same for setting both high and low alarms.

**NOTE:** The alarm setting will not change if the user presses **EXIT** before confirming the change.

## Alarm Conditions


### High Priority

A high priority alarm will alert at 1 second. During a high alarm event, the LED flashes and a high priority alarm tone will sound. Pushing the control knob acknowledges and silences the alarm, but the visual flashing LED continues until the alarm condition resolves. After 30 seconds of a high priority alarm condition, a maximum volume alarm tone will sound, overriding any pre-selected alarm volume level.



### Battery Empty

A **battery empty** high priority alarm will alert at 3% remaining battery power. During this high priority alarm event, the LED flashes and a high priority alarm tone will sound. The alarm can be silenced, but the visual flashing LED continues until the alarm condition is resolved. After 30 seconds of a high priority alarm condition, a maximum volume alarm tone will sound, overriding any pre-selected alarm volume level.



### Loss of Power Alarm

A high priority alarm will sound when the TXP® 5 battery is drained beyond critical, and the device is not connected to AC power. When battery power is less than 2%, the screen shuts off, and the unit emits a high priority alarm tone until the control knob is pressed to silence the alarm tone. The unit will display the battery charge symbol.


**Device will continue to ventilate the patient as set.**



## Medium Priority

### Battery Low

A **low battery** medium priority alarm will alert at 10% remaining battery power. During this medium alarm event, the LED flashes and a medium priority alarm tone will sound. The alarm can be silenced, but the visual flashing LED continues until the alarm condition is resolved. After 30 seconds of a medium priority alarm condition, a maximum volume alarm tone will sound, overriding any pre-selected alarm volume level.



### Patient Pressure Below Set Limit

A medium priority alarm will alert at 1 second. During a medium alarm event, the LED flashes and a medium priority alarm tone will sound. Pushing the control knob acknowledges and silences the alarm, but the visual flashing LED continues until the alarm condition is resolved. After 30 seconds of a medium priority alarm condition, a maximum volume alarm tone overrides any selected volume level.



### Alarm Silence

Pushing the control knob will silence an alarm. Once the alarm is silenced, a 2-minute countdown begins. If, after 2 minutes, the alarm condition persists, normal alarming will resume.

When an alarm is silenced, the alarm screen will flash (on and off) until the alarm condition is resolved.



## Battery Discharge

Battery at 10%	Low Battery Alarm Connect to AC power and begin charging for normal operation.
Battery at 3%	Battery Empty Alarm Connect to AC power and begin charging for normal operation.
Battery at 2%	Loss of Power Alarm The LED is Off. The backlight is off. Battery charge symbol is displayed. The Alarm is sounding. The display is no longer functioning, except the mute button. Once Mute is pressed, it mutes the alarm permanently. Connect to AC power and begin charging for normal operation.
Battery at 1%	Screen is blank; the system is in an extremely low power mode. Connect to AC power to wake the display and begin charging for normal operation.

**NOTE:** Device will operate normally when connected to AC power and charging.


**NOTE:** Ventilation will continue and is unaffected when the display is in low power mode.

## Chapter 5: Setup


### Controller and Stand





### DISS Gas Connection




### Power Supply




## Mounting to Pole/Post




## Adjusting Clamp Range and Tension



## Mounting to Rail



## Phasitron® A50606-TXP Assembly



The Phasitron® A50606-TXP is the mechanical/physiological breathing circuit interface. The Phasitron® has a sliding venturi that acts as both the inhalation and exhalation valve.

1. Connect clear tubing connector (A) to white cap on rear of Phasitron® (B).
2. Connect red tubing connector (C) to front measuring port on Phasitron® (D).
3. Connect Phasitron® delivery port (E) to patient interface.

## Connect Tubing Harness


4. Connect red connector on breathing circuit harness (F) to red bulkhead connector labeled "Monitor."
5. Connect clear connector on breathing circuit harness (G) to silver bulkhead connector labeled "Phasitron."

## Configurations

The Phasitron® A50606-TXP kit can be used invasively or non-invasively using standard endotracheal tubes or mask.


### Intubated Patient with Phasitron®

*Supplemental oxygen can be added after filter.*



### Standard Circuit


Recommended HME/HMEF bacterial/viral filter with an efficiency of  $> 99.999\%$  and 2 cmH<sub>2</sub>O or less flow resistance at 60 L/min.



### Optional PEEP Valve Kit (PRT-A70143)


Optional PEEP valve  
with or without filter

Optional one way valve  
used as balancing  
gate when PEEP valve  
is used



### Optional Secondary Filters

Optional secondary  
filters may be added  
to the two 15 mm OD ports





## Chapter 6: Pre-Use Check


1. Rotate the **Frequency** knob to full clockwise position.



2. Rotate the **Amplitude** knob to full clockwise position.




3. Connect TXP® 5 to gas supply source.



4. Connect the Phasitron® to the TXP® 5 red and silver bulkhead fittings.



5. Close the Phasitron® patient port, with or without filter, ensuring a tight seal.



6. Rotate the **Amplitude** knob to full counterclockwise position.



7. Observe frequency on right side of display less than 250 pulses per minute.

8. Observe Mean Airway Pressure on display greater than 20 cmH<sub>2</sub>O.  
Open the Phasitron® patient port; the disconnect alarm must sound.  
Re-close port to stop alarm and continue.

9. Rotate the **Frequency** knob to the left, full counterclockwise position.



10. Observe an increase in pulse frequency to a rate greater than 550 pulses per minute.

11. Observe Mean Airway Pressure on display greater than 20 cmH<sub>2</sub>O.

12. Rotate the **Frequency** knob to arrow up 12:00 position.



13. Rotate the **Amplitude** knob slowly clockwise. Observe a decrease in amplitude until "off" at the full clockwise position.




14. Check is complete.

## Chapter 7: Ventilation Settings

### Standard Circuit Initial Settings

1.	Set Frequency by turning arrow straight up (rate of 500 +/- 25).
2.	Attach Phasitron® to patient hospital approved HMEF.
3.	Start with Amplitude control knob turned full clockwise (right).
4.	Slowly turn Amplitude control knob counterclockwise (left), until patient's chest is observed to be moving (wiggle). Ensure that the chest is moving, just below the ribs.
5.	Observe patient SpO <sub>2</sub> and CO <sub>2</sub> .
6.	After 30 minutes, draw an ABG.
7.	Record Mean Airway Pressure (MAP), high-frequency rate, and amplitude along with SpO <sub>2</sub> and CO <sub>2</sub> .

### Adjustment Options: Standard Circuit



#### Increase pO<sub>2</sub> and decrease CO<sub>2</sub>

- Increase amplitude in steps of 2-4 cmH<sub>2</sub>O.


#### Increase O<sub>2</sub> if CO<sub>2</sub> is OK

- If device is plugged into an O<sub>2</sub> outlet, the delivered FiO<sub>2</sub> will be 60%.
- If more than 60% FiO<sub>2</sub> is needed, oxygen may be added between the Phasitron® and the patient to achieve close to 100%.
- Increase Frequency by 100 (repeat to a maximum of 700). This may increase MAP which may affect pO<sub>2</sub>.

#### Decrease CO<sub>2</sub> if O<sub>2</sub> is OK

- Increase Amplitude in steps of 2-4 cmH<sub>2</sub>O (while keeping MAP within desired range).
- Decrease Frequency by 100 pulses in a stepwise fashion to a low of 300.

## Adjustment Options: Standard Circuit with Optional PEEP Valve and Inspiratory Valve



### Increase $pO_2$ and Decrease $CO_2$

- Increase Amplitude in steps of 2-4  $cmH_2O$ . Watch for increase in MAP.

### Increase $O_2$ if $CO_2$ is OK


- If the device is plugged into an  $O_2$  outlet, the delivered  $FiO_2$  will be 60%.
- If more than 60%  $FiO_2$  is needed, oxygen may be added between the Phasitron® and the patient to achieve close to 100%.
- Increase the Mean Airway Pressure (MAP) by 2  $cmH_2O$  with the mechanical PEEP valve.
- Increase Frequency by 100 (repeat to a maximum of 700); this may increase  $CO_2$ .

### Decrease $CO_2$ if $O_2$ is OK

- Increase Amplitude in steps of 2-4  $cmH_2O$  (while keeping MAP within a desired range).
- Decrease Frequency by 100 pulses, in a stepwise fashion, to a low of 300.

### Increase $CO_2$ if $O_2$ is OK


- Decrease Amplitude in steps of 2-4  $cmH_2O$  (while keeping MAP within a desired range). If MAP is decreasing, adjust PEEP valve to a higher PEEP setting for similar MAP.
- Increase Frequency by 100 pulses to a maximum of 700.



Amplitude


High-frequency rate


MAP



## Patient Monitoring

Clinicians should evaluate how their patients tolerate the ventilation. Auscultation and observation of the mechanical vibrations of the chest and abdomen are primary indicators of effective treatment. During HFPV, it is vital to maintain an unobstructed and unrestricted airway.

 **WARNING:** Follow proper suctioning procedures to maintain a patent airway.

 **WARNING:** Regular patient assessment, along with continuous monitoring of SpO<sub>2</sub> and end-tidal CO<sub>2</sub>, are necessary to ensure blood gases are at the proper level.

## Assessment

• SpO <sub>2</sub>	• Blood pressure and perfusion
• Chest wiggle	• Breath sounds
• Ventilator settings	• Phasitron® and ETT positioning
• Blood gas ~ 30 minutes after starting HFPV or after an adjustment and then PRN per institution protocol	

## Ventilation

Control over arterial pCO<sub>2</sub> is achieved by manipulating delivered volumes.

- $HF = VCO_2 + (VT)2 \times f$

## Increase Ventilation

To increase ventilation and decrease pCO<sub>2</sub>:

- Increase amplitude in 2-4 cmH<sub>2</sub>O increments. Amplitude is directly related to volume delivery.
- Decrease frequency. This will increase pulse to pulse time, therefore increase volume.

## Oxygenation

Control over arterial  $pO_2$  is achieved by manipulating  $FiO_2$  or increasing the Mean Airway Pressure. The following will increase  $FiO_2$  and/or Mean Airway Pressure, increasing arterial  $pO_2$ .

### Increase Oxygenation

<b>Increase <math>FiO_2</math></b>	Connect the unit to a blended air/ $O_2$ system; adjust $FiO_2$ on blender. If the TXP® 5 is connected to an oxygen source between the patient and the filter, $FiO_2$ of approximately 98% can be expected. See page 31 for $FiO_2$ measurements.
<b>Increase Amplitude</b>	This will increase Mean Airway Pressure, which will increase $pO_2$ . There will also be an increase in volume, which will decrease $pCO_2$ .
<b>Increase Frequency</b>	This will increase Mean Airway Pressure, which will increase $pO_2$ . This may also increase $pCO_2$ .
<b>Attach accessory PEEP valve to Exhalation port and one-way balancing gate valve to the Inhalation port</b>	This will increase FRC, which will increase Mean Airway Pressure. This will increase MAP, without changing Amplitude.

In the face of a challenging clinical scenario with worsening ventilation and/or oxygenation, a complete assessment of the patient is necessary to find the cause for V/Q inequality. Common V/Q problems associated with acute or chronic respirator injuries are noted below.


**NOTE:** Any increase in MAP will increase  $pO_2$ . Any increase in Amplitude or Frequency will increase MAP.

**NOTE:** Amplitude and Frequency are interconnected, changing one parameter will affect the other.

### Decrease Ventilation: Increase $pCO_2$

- Reduce amplitude (this will also decrease MAP): Decrease in 2-4 cmH<sub>2</sub>O increments.
- Frequency: Increase in increments of 60.
- Mechanical deadspace: Add 6-inch corrugated tubing (as shown below).

**Mechanical deadspace**




**NOTE:** In the face of a challenging clinical scenario with worsening ventilation and/or oxygenation, complete assessment of the patient to find the cause of V/Q inequality.

## Chapter 8: Cleaning and Maintenance

### Cleaning

**NOTE:** All single-patient use components and Phasitron® are not intended for cleaning, sterilization, or re-use. Replace single-patient use components regularly, following your healthcare institution's protocol.

Before cleaning any part of the TXP® 5, disconnect external power sources.

 **WARNING:** Do not perform maintenance or service on the TXP® 5 while it is powered on or in use. Maintenance or service procedures performed during use may temporarily alter the performance and result in patient harm.


### TXP® 5 Controller

Clean the controller according to hospital/institutional protocols. Always clean between patients and when visibly soiled. Clean the controller with a clean, lint-free cloth or paper towel moistened with water (including water mixed with soap or a mild detergent) or 70% isopropyl alcohol.

 **CAUTION:** Do not spray any cleaning solution directly onto the controller.


 **CAUTION:** Do not immerse or allow liquids to access the controller.


 **CAUTION:** Use only approved cleaning solvents recommended in this manual.

 **CAUTION:** Before plugging it in to an external power source, allow the TXP® 5 to dry completely after cleaning.

### Digital Display Screen

Clean the display screen using 70% isopropyl alcohol and according to facility protocols. Do not spray any type of cleaner directly onto the LCD.

 **CAUTION:** The use of cleaning methods not outlined in these instructions may cause damage to the display.

 **WARNING:** The cell used in this device may present a risk of fire or chemical burn hazard if mistreated. Do not disassemble, heat above 100°C (212°F), or incinerate. The use of another cell may present a risk of fire or explosion.





*At the end of useful life of a TXP® 5 unit, disposal should be in accordance with local, state, federal, and international laws.*

## Phasitron® A50606-TXP Breathing Circuit Kit

The Phasitron® A50606-TXP is a single patient device. Follow hospital guidelines for cleaning and disinfection. Percussionaire® recommends changing the Phasitron® breathing circuit every seven days or sooner if visibly soiled. Do not exceed seven days.

 **WARNING:** The Phasitron® kit is for single patient use only.

 **CAUTION:** Do not use harsh cleaners, solvents or detergents. Equipment damage could occur.

 **CAUTION:** Failure to follow the manufacturer's cleaning instructions could cause equipment damage.

## Cleaning and Disinfecting Solutions

The Phasitron® breathing circuit has been tested for biocompatibility with the following cleaning and disinfecting solution:

Chemical Class	Active Ingredient
Bleach	5.25% Sodium hypochlorite
Alcohol	70% Isopropyl alcohol
Peroxide	3% Hydrogen peroxide chloride

## Maintenance

### TXP® 5 Controller


The TXP® 5 should be checked by an authorized Percussionaire® service technician anytime clinical efficacy is not as expected. Perform a pre-use check before ventilating a patient. Perform preventive maintenance and functional evaluation annually.

### Lithium-ion Battery Care

The internal battery is lithium-ion. To maximize battery life, charge the battery before it drains completely. The battery will charge whenever the TXP® 5 is plugged into an external source of power.


Over time, the capacity of the battery will diminish. Percussionaire® recommends that you replace the battery when it no longer holds a charge. The internal battery must be replaced by authorized Percussionaire® service technicians or Percussionaire®-certified biomed technicians.

## Battery Replacement




The TXP® 5 uses a Lithium-ion 3.7V 2600mAh rechargeable battery (18650). Percussionaire® part #B21111

**NOTE:** Replace the battery every 3 years .



Open the top case by removing 4 screws.

Use a small flat screwdriver to remove the safety clip from battery.



Insert a new battery noting the positive terminal. Replace the safety clip.

**NOTE:** Ensure the battery is fully charged before storing the TXP® 5.

**NOTE:** Only use batteries provided by Percussionaire®.

**Do not use batteries from other manufacturers.**

**⚠ WARNING:** Do not damage the rechargeable lithium-ion battery. A damaged battery may cause an explosion or fire and may result in personal injury and/or property damage. To prevent injury or damage:


- Do not use or charge the battery if it appears to be damaged. Signs of damage include, but are not limited to, discoloration, warping, and leaking battery fluid.
- Do not expose the battery to fire, high temperature.
- Do not immerse the battery in water.
- Do not use or store the battery inside a vehicle during hot weather.
- Do not drop or puncture the battery.
- Do not open the battery or short-circuit its contacts.

**⚠ WARNING:** Avoid contact with the rechargeable lithium-ion battery if it appears to be leaking. Battery fluid is corrosive. Contact with battery fluid may result in personal injury and/or property damage. To prevent injury or damage:

- If the battery leaks, avoid contact with the battery fluid.
- If the battery fluid gets into your eyes, immediately rinse your eyes with clean water and seek medical attention. Do not rub your eyes.
- If battery fluid gets onto your skin or clothing, immediately use clean water to wash off the battery fluid.




## Chapter 9: Troubleshooting

 **WARNING:** If there are any unexplained changes in the performance of the device, if the device makes unusual sounds or is damaged in any way, discontinue use. Begin troubleshooting process and contact distributor or an authorized Percussionaire® service center.

Problem	Possible Cause	Corrective Action
No Percussion	No air or gas supply to the TXP® 5	Confirm gas supply can deliver 50 psig pressure at 15 LPM flow rate. Connect to a different gas source to confirm pressure and capacity.
	TXP® 5 has air or gas supply, but no percussion.	Make sure the Amplitude knob is turned fully counterclockwise, and the Phasitron® red and clear tubing is connected
	Internal percussion valve failed.	Send device to an authorized service center for repair.
Low Amplitude	Low air or gas supply pressure/flow capacity	Confirm gas supply connection can deliver 50 psig pressure at 15 LPM flow rate. Connect to a different gas source to confirm pressure and capacity.
	Gas supply air leak	Listen for air leaks. Visually inspect both the gas supply connections and the Phasitron® to TXP® 5 connections. Inspect for damaged O-rings at the tubing connectors.
	Correct amplitude, but the display is indicating a lower than expected pressure.	The TXP® 5 displays the pulse amplitude pressure in the small upper-right part of the display. The numbers near the center of the display indicate Mean Airway Pressure. This is a lower number than the pulse amplitude pressure. Check with a test lung.

Continued on page 30

Problem	Possible Cause	Corrective Action
Frequency does not change	The Frequency knob is loose on control shaft.	<p>The Frequency knob must be turned fully counterclockwise, to the stop. Remove the Frequency knob with a 3/32" hex key wrench and rotate the valve shaft for a frequency of 200 ppm. Reinstall the knob and confirm the frequency is about 200 ppm at the full clockwise position.</p> 
TXP® 5 has no display	No signal to display	The Phasitron® must be connected to the clear and red tubing. Silver and red connectors on the TXP® 5 face panel must be connected to the Phasitron®.
	Display not activated	The Phasitron® patient port must see a pressure of 2.5 cmH <sub>2</sub> O for more than one second to enter wake mode.
	Display does not stay on	Confirm the Phasitron® is connected to the patient, test lung, or occluded for more than 5 seconds to activate display.
	Display not functional	Check power connection and reconnect
TXP® 5 has wrong display	Power-On Self-Test and calibration were unable to start correctly	Turn off device and unplug power. Wait 5 seconds and reconnect power.
	Failed display software	Send the display to an authorized service center for repairs.
Device performance changes	Electrical interference	Move TXP® 5 away from any potential sources of electromagnetic interference (EMI) including MRI equipment, medical imaging systems, security systems, appliances, wireless communications equipment (such as cellular phones), computers, and televisions.

## Chapter 10: Technical Specifications

### TXP® 5 Ventilator

<b>Accessories</b>	Phasitron® Kit A50606-TXP
<b>Pulse/Interval Ratio</b>	Frequency: 250-350    Ratio: 1:2 Frequency: 350-450    Ratio: 1:1.5 Frequency: 450+        Ratio: 1:1
<b>Operating Range</b>	Temp., 0°C to 49°C (32°F to 120°F) Humidity 5%-95%
<b>Storage and Transport Range</b>	Temp., -20°C to 60°C (-4°F to 140°F) Humidity < 93% non-condensing
<b>Power Input</b>	90-264 VAC, 50/60Hz, 1 amp
<b>Battery Type</b>	Li-ion 18650 3.7V 2600mAh rechargeable
<b>Run Time</b>	Continuous
<b>Pulse Frequency</b>	200-700 pulses per minute (approximately)
<b>High/Low Pressure Settings</b>	Digital display
<b>Mean Airway Pressure</b>	Digital display, 0-99 cmH <sub>2</sub> O/hPa
<b>Pulse Amplitude</b>	Digital display, 0-99 cmH <sub>2</sub> O/hPa Accurate to +/- 1 cmH <sub>2</sub> O
<b>Alarm: Airway High-Pressure Limit</b>	Audible indicator Red LED display during audible alarm
<b>Alarm: Airway Low-Pressure Limit</b>	Audible indicator Amber LED display during audible alarm
<b>Alarm: Circuit Disconnect Alarm</b>	Audible indicator Red LED display during audible alarm
<b>Gas Source: Hospital Wall Gas, Compressed Oxygen Tanks, or Compressor Flow</b>	50-80 psi, 3.45-5.5 BAR
<b>Flow</b>	15 LPM average
<b>Dimensions (W x H x D)</b>	15.875 cm H x 11.43 cm W x 15.24 cm D (6.25 in H x 4.5 in W x 6.0 in D)
<b>Weight</b>	.5 kg (1.16 lbs)
<b>Required Maintenance</b>	Annual function check/Replace the battery every 3 years

### Measured Performance

TXP <sup>5</sup> Frequency		R=50, C=10 mL/cmH <sub>2</sub> O		R=20, C=10 mL/cmH <sub>2</sub> O		R=5, C=10 mL/cmH <sub>2</sub> O	
300	MAP Amplitude	40 cmH <sub>2</sub> O 99+cmH <sub>2</sub> O	31 LPM	30 cmH <sub>2</sub> O 79 cmH <sub>2</sub> O	56 LPM	21 cmH <sub>2</sub> O 47 cmH <sub>2</sub> O	78 LPM
500	MAP Amplitude	47 cmH <sub>2</sub> O 99+cmH <sub>2</sub> O	28 LPM	39 cmH <sub>2</sub> O 72 cmH <sub>2</sub> O	48 LPM	32 cmH <sub>2</sub> O 54 cmH <sub>2</sub> O	67 LPM
700	MAP Amplitude	52 cmH <sub>2</sub> O 99+cmH <sub>2</sub> O	23 LPM	48 cmH <sub>2</sub> O 87 cmH <sub>2</sub> O	36 LPM	44 cmH <sub>2</sub> O 66 cmH <sub>2</sub> O	50 LPM

Sampled from Phasitron® breathing circuit PN A50606-TXP, inline with TSI 4000 Flow Meter, inline with Ingmar Medical's QuickLung®.

<sup>1</sup>Resistance on QuickLung® set at each of the available resistance settings, while compliance remained set at 10 ml/cmH<sub>2</sub>O.

<sup>2</sup>Frequency recorded from TXP® 5 display.

FiO <sub>2</sub> Measurements					
<b>TXP<sup>5</sup></b> Frequency		R=50, C= 10 mL/cmH <sub>2</sub> O		O <sub>2</sub> Sat without added O <sub>2</sub>	O <sub>2</sub> Sat with added O <sub>2</sub> post filter
300	MAP/AMPLITUDE	32/99+	31 LPM	67	94.6 @ 5 LPM
500	MAP/AMPLITUDE	37/92	28 LPM	67	100 @ 5 LPM
700	MAP/AMPLITUDE	39/89	25 LPM	66	100 @ 5 LPM
<b>TXP<sup>5</sup></b> Frequency		R=20, C= 10 mL/cmH <sub>2</sub> O		O <sub>2</sub> Sat without added O <sub>2</sub>	O <sub>2</sub> Sat with added O <sub>2</sub> post filter
300	MAP/AMPLITUDE	23/69	54 LPM	59.7	95 @ 8.5 LPM
500	MAP/AMPLITUDE	29/63	50 LPM	66.5	96.1 @ 8 LPM
700	MAP/AMPLITUDE	34/69	42 LPM	72.1	95.1 @ 4.5 LPM
<b>TXP<sup>5</sup></b> Frequency		R=5, C= 10 mL/cmH <sub>2</sub> O		O <sub>2</sub> Sat without added O <sub>2</sub>	O <sub>2</sub> Sat with added O <sub>2</sub> post filter
300	MAP/AMPLITUDE	14/33	80 LPM	55.8	91.5 @15 LPM
500	MAP/AMPLITUDE	20/38	75 LPM	62.3	100 @ 10 LPM
700	MAP/AMPLITUDE	25/46	68 LPM	73	96.1 @ 8 LPM

Sampled from Phasitron® breathing circuit PN A50606-TXP, inline with TSI 4000 Flow Meter, inline with Ingmar Medical's QuickLung®.

<sup>1</sup>Resistance on QuickLung® set at each of the available resistance settings, while compliance remained set at 10ml/cmH<sub>2</sub>O.

<sup>2</sup>Frequency recorded from TXP® 5 display.

## Digital Display Specifications

<b>Storage and Transport Range</b>	Temp., -20°C to 60°C (-4°F to 140°F) Humidity <93% non-condensing
<b>Operating Range</b>	Temp., -20°C to 60°C (-4°F to 140°F) Humidity <93% non-condensing
<b>Display:</b>	
Display Type	Transreflective
Visible Area	58.00 mm W x 28.80 mm H
Backlight	LED-white
Pixels	128 x 64
Graphics Color	Black (White – inverted)
Background Color	White (Black – inverted)
Spot Size	0.36 mm W x 0.36 mm H
Step of Points	0.40 mm H x 0.40 mm H
<b>Rate Range</b>	200-700 pulses per minute (approximately)
<b>Pressure Range</b>	0-99 cmH <sub>2</sub> O/hPa
<b>Pressure Accuracy</b>	Accurate to +/- 1 cmH <sub>2</sub> O
<b>Battery</b>	Percussionaire® part #B21111, Lithium-ion 18650 3.7V

## Phasitron® A50606-TXP

<b>Size</b>	13.5 mm x 17 mm (5 ¼" x 6 ¾")
<b>Weight</b>	123 g (0.27 lb)
<b>Operating Range</b>	Temp., 0°C to 49°C (32°F to 120°F) Relative humidity range 5% to 95%
<b>Storage and Transport Range</b>	Temp., -20°C to 60°C (-40°F to 140°F) < 93% non-condensing
<b>Rate Range</b>	0-999 pulses per minute
<b>Pressure Range</b>	0-150 cmH <sub>2</sub> O/hPa
<b>Filtration Efficiency</b>	BFE 99.999%, VFE 99.9999%
<b>Red Line Filter</b>	0.027 micron hydrophobic
<b>Disposal</b>	Recycle according to local laws
<b>Service Life</b>	7 days
<b>Shelf Life</b>	2 years from date of manufacture

## Battery: Percussionaire® part # B21111

<b>Battery Type</b>	Lithium-ion 18650 3.7V 2600mAh rechargeable
<b>Diameter</b>	19 mm
<b>Length</b>	69.7 mm
<b>Rated Capacity Nominal</b>	2600mAh
<b>Discharge Cut-Off Voltage</b>	2.75V
<b>Charging Voltage</b>	4.2V
<b>Shipments Voltage</b>	3.5-3.65V SOC <= 30%
<b>Impedance</b>	≤100 mΩ
<b>Standard Charge:</b>	
<b>Constant Current</b>	.2C
<b>Constant Voltage</b>	4.2V
<b>Cut-off Current</b>	0.02C
<b>Standard Discharge:</b>	
<b>Constant Current</b>	.2C
<b>End Voltage</b>	2.75V
<b>Maximum Charge Current:</b>	
<b>Constant Current</b>	0.5C
<b>Constant Voltage</b>	4.2V
<b>Cut-off Current</b>	0.02C
<b>Fast Discharge:</b>	
<b>Constant Current</b>	0.5C
<b>End Voltage</b>	2.75V
<b>Maximum Continuous Discharge Current</b>	3.5A
<b>Operation Charge</b>	0 to 45 Celsius, 60 +/- 25% R.H.
<b>Temperature Range Discharge</b>	Temperature Range Discharge
<b>Cycle life</b>	300+ Cycles
<b>Storage Temperature:</b>	
<b>&lt;= 1 month</b>	-20 to 50 Celsius
<b>&lt;= 3 months</b>	-10 to 30 Celsius
<b>&lt;= 1 year</b>	0 to 30 Celsius
<b>Weight</b>	50 g (approximate)

**NOTE:** Only use batteries provided by Percussionaire®.  
Do not use batteries from other manufacturers.

## Chapter 11: Service and Repair

Percussionaire® Corporation recommends an annual preventive maintenance (PM) for the TXP® 5. An annual PM consists of a thorough cleaning, functional evaluation, and if necessary, a recalibration.

The standard interface connector (USB) is provided on the internal PCB, which is non-active and only used for calibration and software update. This device does not allow communication.

To return a device to a factory service center for repair, contact [customerservice@percussionaire.com](mailto:customerservice@percussionaire.com) or call (208) 263-2549.

## Disposal of Equipment



At the end of useful life of a TXP® 5 unit, disposal should be in accordance with local, state, federal, and international laws.

## Battery Disposal

The battery is considered electronic waste and must be disposed of according to local regulations. Follow local governing ordinances and recycling plans regarding disposal or recycling of the battery.



At the end of useful life of a TXP® 5 unit, disposal should be in accordance with local, state, federal, and international laws.

## Chapter 12: Limited Warranty

Percussionaire® warrants that the TXP® 5 shall be free from defects of workmanship and materials and will perform per the product specifications for five years from the date of purchase (proof of delivery will be required). The device must be run for 30 minutes annually if kept in storage conditions. Batteries carry a 12-month warranty. If the product fails to perform in accordance with the product specifications, Percussionaire® will repair or replace – at its option – the defective material or part. Percussionaire® will pay customary freight charges to and from Percussionaire® or an authorized Percussionaire® service center. This warranty does not cover damage caused by non-approved cleaning or sterilization, accident, misuse, abuse, alteration, or other defects not related to material or workmanship. Percussionaire® disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product.



**PERCUSSIONNAIRE®**  
**CORPORATION**

🏠 130 McGhee Road., Suite 109, Sandpoint, Idaho 83864 USA

📧 [percussionnaire.com](http://percussionnaire.com) ☎ +1.208.263.2549