EN   Cleaning and Disinfection .....	2
DA   Rengøring og desinfektion .....	7
DE   Reinigung und Desinfektion.....	12
ES   Limpieza y desinfección .....	17
FI   Puhdistus ja desinfiointi .....	22
FR   Nettoyage et désinfection .....	27
IT   Pulizia e disinfezione .....	32
JA   清掃および消毒 .....	37
NL   Reiniging en desinfectie.....	42
NO   Rengjøring og desinfeksjon .....	47
PL   Czyszczenie i dezynfekcja .....	52
PT   Limpeza/Desinfeção.....	57
SV   Rengöring och desinficering.....	62

## EN | Cleaning and Disinfection

The cleaning and disinfection solutions listed here have been tested for their suitability for the following products:

Product Name	Referred to as	REF
Sentec Digital Monitor	Monitor	SDM
Digital Sensor Adapter Cable	Cable	AC-150, AC-250, AC-750
tCOM+ Transcutaneous Monitor	Monitor	103164
Sensor Adapter Cable	Cable	103420, 103421, 103422
V-Sign™ Sensor 2	Sensor (Including sensor and cable)	VS-A/P/N
OxiVenT™ Sensor	Sensor (Including sensor and cable)	OV-A/P/N

The following list of name brand products for cleaning and disinfection may be substituted by other name brand products of equivalent composition. However, as there are many factors (contact time, force applied during wiping etc.) that can influence the functionality of the sensors, cables, and monitors, it is the responsibility of the user to test the functionality of the sensors, cables, and monitors after cleaning and/or disinfection.

**NOTE:** For the specific use, handling, storage, and other recommendations specific to the disinfection solutions listed here, follow the manufacturer's instructions.

**Sentec recommends using 70% Isopropanol for surface cleaning and low-level disinfection, and as a final wipe to remove residues after disinfection with the other solutions.**

Table 1: Surface Cleaning and Low-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
WEBCOL™ Alcohol Prep Pads (Cardinal Health™/Covidien™)	70% Isopropanol.	Carefully wipe all the surfaces of sensors, cables, and monitors. Keep damp for 1 minute. Allow to dry.
mikrozid® universal wipes/universal wipes premium (Schülke&Mayr)	17.4 g Propan-2-ol; 12.6 g Ethanol (94% w/w); <5% anionic surfactants.	
mikrozid® AF wipes (Schülke&Mayr)	25 g Ethanol (94%); 35 g Propan-1-ol.	

Table 2: Low-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
<b>kodan® forte</b> (Schülke&Mayr)	45.0% 2-Propanol; 10.0% 1-Propanol; 0.20% Biphenyl-2-ol; Hydrogen peroxide solution (30%).	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes soaked in Kodan® forte. Wait 5 min. Remove all Kodan® forte residues using 70% Isopropanol. Allow to dry.
<b>gigasept® AF</b> (Schülke&Mayr)	15 g Didecylmethylammonium chloride; 15 g Phenoxypropanole; 6.9 g Aminoalkyl glycines; 15-30% non-ionic surfactants.	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes moistened with 4% solution (e.g., 960 ml deionized water plus 40 ml gigasept® AF). Wait 15 minutes. Remove all gigasept® residues using 70% Isopropanol. Allow to dry.
<b>terralin® protect</b> (Schülke&Mayr)	22 g Quaternary ammonium compounds; Benzyl-C12-16-alkyldimethyl chloride; 17 g 2-phenoxyethanol; 0.9 g amines; C12-14 alkyldiaminoethylglycine HCl; 5 - 15% non-ionic surfactants.	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes moistened with 2% solution (e.g., 980 ml deionized water plus 20 ml terralin® protect). Wait 15 minutes. Remove all terralin® residues using 70% Isopropanol. Allow to dry.  <b>NOTE:</b> terralin® protect can cause micro cracks if used regularly, over time on a SDM or tCOM+ Transcutaneous Monitor and damage sensors and cables. Always remove all residues from surfaces using 70% Isopropanol.
<b>Sanoclean AR</b> (Sanosil)	1.5% Hydrogen peroxide (0.003% silver).	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes soaked in Sanoclean AR. Wait 15 minutes. Remove all residues using 70% Isopropanol. Allow to dry.  <b>NOTE:</b> The commercial solution is ready to use, no further dilution required.
<b>Dismozon® plus</b> (Bode Chemie)	95.8% Magnesium monoperoxyphthalate hexahydrate.	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes moistened with 3.6% solution (e.g., 36 g of Dismozon® granulate per liter of deionized water). Wait 15 minutes. Remove all Dismozon® residues using 70% Isopropanol. Allow to dry.
<b>Mikrobac® Forte</b> (Bode Chemie)	20% Benzyl-C12-18-alkyldimethylammoniumchloride 5% N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine.	Carefully wipe all the surfaces of sensors, cables, and monitors with wipes moistened with 2% solution (e.g., 980 ml deionized water plus 20 ml Mikrobac® Forte). Wait 15 minutes. Remove all Mikrobac® Forte residues using 70% Isopropanol. Allow to dry.  <b>NOTE:</b> Mikrobac® Forte can cause micro cracks if used regularly, over time on a SDM or tCOM+ Transcutaneous Monitor and damage sensors and cables. Always remove all residues from surfaces using 70% Isopropanol.

Table 2: Low-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
Super Sani-Cloth® (PDI)	0.25% n-Alkyldimethylethylbenzyl-ammonium chlorides; 0.25% n-alkyl dimethyl benzyl ammonium chlorides; 55% Isopropyl Alcohol.	Carefully wipe all the surfaces of sensors, cables, and monitors with the pre-moistened wipe. Wait 2 minutes. Remove all Super Sani-Cloth® residues using 70% Isopropanol. Allow to dry.   <b>CAUTION:</b> This cleaning agent is potentially aggressive to the skin. Thus, change the sensor membrane after cleaning to protect the skin.
CaviCide™ or CaviWipe™ (Metrex® research corporation)	<20% Isopropanol; <5% 2-Butoxyethanol; < 1% Benzo/Didecyldimethyl Ammoniumchlorid.	Spray CaviCide™ or use the CaviWipe™ directly on the surfaces of sensors, cables, and monitors.  Allow surface to remain visibly wet for 3 minutes. Remove all CaviCide™ or CaviWipe™ residues using 70% Isopropanol. Allow to dry.   <b>CAUTION:</b> This cleaning agent is potentially aggressive to the skin. Thus, change the sensor membrane after cleaning to protect the skin.
Prevantics® Device Swab (PDI)	3.15% Chlorhexidine Gluconate (w/v); 70% Isopropyl Alcohol (v/v).	Carefully wipe all the outer surfaces with the pre-moistened PDI wipe for 5 seconds. Allow to dry for 5 seconds. Remove all PDI residues using 70% Isopropanol. Allow to dry.

Table 3: High-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
Sani-Cloth® Bleach (PDI)	Sodium hypochlorite: Approx. 6000 ppm / 5500ppm available free chlorine.	Carefully wipe all the surfaces of sensors, cables, and monitors. Wait 4 minutes. Remove all residues using 70% Isopropanol. Allow to dry.   <b>CAUTION:</b> Do not use bleach cleaners on sensors without a membrane or with a defective membrane. This may damage the PCO <sub>2</sub> unit.  <b>CAUTION:</b> Bleach is potentially aggressive to the skin. Thus, change the sensor membrane after cleaning to protect the skin.
Clorox Healthcare® Bleach Germicidal Wipes	0.55% Sodium hypochlorite; 0.52% (5200 ppm) available chlorine.	 <b>CAUTION:</b> Bleach can corrode metal. Therefore, use bleach cleaners on outer surfaces only and do not bring into contact with metallic parts. Always remove all residues using 70% Isopropanol.

Table 3: High-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
Cidex OPA® (J&J)  <b>NOTE:</b> Requires pre-cleaning using an enzymatic cleaner: ENZOL® (J&J) or Prolystica® (Steris); Subtilisin	0.55% ortho-Phthalaldehyde; Subtilisin.	Sensor procedure: <ol style="list-style-type: none"> <li>1. Remove debris using an enzymatic detergent: Soak the sensor and the desired length of cable in enzymatic solution for 5 minutes. Rinse the sensor and the cable briefly in deionized water. Blot dry carefully.</li> <li>2. Soak the sensor and the desired length of cable in CIDEX OPA solution for 12 minutes.</li> <li>3. Soak the sensor and the desired length of cable 3 times for 2 minutes in sterile or deionized water (replace after each soaking cycle).</li> <li>4. Carefully wipe all surfaces of the sensor and cable with 70% Isopropanol. Allow to dry.</li> </ol> <p><b>!</b> <b>CAUTION:</b> Change the sensor membrane after cleaning to ensure optimal sensor monitoring performance</p> Monitor procedure: <ol style="list-style-type: none"> <li>1. Remove debris using an enzymatic detergent: Wipe the outer surfaces with a wipe moistened with enzymatic solution and wait 5 minutes. Remove enzymatic cleaner by wiping all surfaces with towel moistened in water.</li> <li>2. Perform high level disinfection: Wipe the outer surfaces with a wipe moistened in Cidex OPA® and wait 12 minutes. Remove Cidex OPA® solution by wiping all surfaces 3 times with towel moistened in water. Use a fresh wipe each time.</li> <li>3. Final wipe: Remove all residues using 70% Isopropanol. Allow to dry.</li> </ol>
mikrozyd® PAA wipes (Schülke&Mayr)	0.06% Per-Acetic acid; Hydrogen Peroxide; Acetic acid.	Carefully wipe all the surfaces of sensors, cables, and monitors. Wait 15 minutes. Remove all Mikrozyd PAA residues using 70% Isopropanol. Allow to dry.

Table 3: High-level Disinfection		
Name Brand / Manufacturer	Composition	Recommended Procedure
The following agents are only suitable for sensors and cables:		
Bacillol® AF (Bode Chemie GmbH)	45% propan-1-ol; 25% propan-2-ol; 4.7% ethanol.	Carefully wipe all the surfaces of sensors and cables with wipes soaked in Bacillol® AF  <b>NOTE:</b> The commercial solution is ready to use, no further dilution required.  Allow to dry. Remove all Bacillol® AF residues using 70% Isopropanol. Allow to dry.
Incidin™ Pro (Ecolab Healthcare)	<0.5% Benzalkonium chloride; <0.1% 2-phenoxyethanol; <0.08% N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine.	Carefully wipe all the surfaces of sensors and cables with wipes soaked in 1% solution (e.g. 990 ml deionized water plus 10 ml of Incidin™ Pro). Wait 15 minutes.  Remove all Incidin™ Pro residues using 70% Isopropanol. Allow to dry.  <b>NOTE:</b> Disinfectants containing amine or amine derivatives shall not be used interchangeably with disinfectants containing aldehyde. Thoroughly clean when changing disinfectant.
Clinell® Universal Wipes (GAMA Healthcare)	≤0.5% Benzalkonium chloride; ≤0.5% Didecylmethyl ammonium chloride; ≤0.10% Polyhexamethylene biguanide /PHMB).	Carefully wipe all sensors and cables with the presoaked wipe for at least 1 minute. Leave to air dry.  Remove all Clinell® residues using 70% Isopropanol. Allow to dry. It is not mandatory to change the membrane after cleaning.

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