

Continuous, noninvasive tcPCO₂ monitoring of patients with chronic respiratory failure using the SenTec Digital Monitoring System

*Dr. Friederike S. Magnet, Department of Pneumology,
Cologne-Merheim Hospital, Kliniken der Stadt Köln gGmbH*

Introduction

The pulmonary system comprises two different aspects: the lung and the respiratory pump [1]. A failure of the lung results in hypoxemic respiratory failure (type 1 respiratory failure), whereas a failure of the respiratory pump results in hypercapnic respiratory failure (type 2 respiratory failure) [1]. PCO₂ measurement has become part of the standard of care in the diagnosis and treatment of respiratory failure [1]. Pulse oximetry alone is unable to adequately detect hypoventilation, particularly in patients with healthy lungs, where no ventilation perfusion mismatch occurs [2-5]. In addition to day-time measurements, it is crucial that patients with borderline hypercapnia are also monitored during the night, as nocturnal hypoventilation might be missed if only daytime partial pressure of carbon dioxide (PCO₂) measurements are performed [6].

There are invasive and noninvasive methods for monitoring PCO₂. Blood gas analysis is the gold standard and provides more than just PCO₂ values, but it has the drawback of being an invasive technique, which could cause the patient pain and could therefore falsify sleep studies [1]. It is also only able to provide a spot-check analysis. By contrast, transcutaneous PCO₂ monitoring (tcPCO₂) allows continuous noninvasive PCO₂ measurement, is painless and does not disrupt the patient's sleep at night [7].

Clinical evidence

Numerous trials have evaluated the use of tcPCO₂ in different clinical applications.

Noninvasive ventilation

tcPCO₂ has proven to be reliable in predicting PCO₂ during the initiation of acute noninvasive ventilation (NIV) therapy in COPD patients as well as in nocturnal monitoring of noninvasive ventilated patients with type 2 chronic respiratory failure [7-9].

It was shown that tcPCO₂ measurement is able to reduce blood gas analysis during NIV initiation and that tcPCO₂ values measured two minutes after a blood gas analysis were the best predictor of PaCO₂ (=lag time) [9]. Furthermore, tcPCO₂ was able to detect nocturnal hypoventilation during NIV check-ups more effectively than blood gas analysis, which only provides a spot-check of the respiratory situation (see figure 1). Continuous overnight tcPCO₂ monitoring (SenTec Monitor System) revealed variations with a median of 12.3 mmHg, whereas capillary PaCO₂ measurements showed variations with a median of only 6.3 mmHg [7]. German guidelines recommend tcPCO₂ monitoring combined with blood gas analysis during the initiation of NIV [1, 6].




Figure 1: Overnight trend of PaCO₂ (black boxes) and tcPCO₂ (orange line) showing nighttime fluctuations of PCO₂ [7]

Acute respiratory failure

tcPCO₂ was also evaluated in a mixed cohort of patients with acute respiratory failure (ARF). It was shown to be superior to end-tidal CO₂ (etCO₂) in predicting PaCO₂ in non-intubated ICU patients (primarily with acute exacerbated COPD, but also ARDS, heart failure etc.) [10]. Another trial used tcPCO₂ in patients with ARF in the emergency department and demonstrated a significant correlation to PaCO₂ ($R^2=0.83$, $p<.001$), with better agreement in normothermia compared to hyperthermia as well as PaCO₂ values < 56 mmHg compared to > 60 mmHg [11].

Nasal high-flow oxygen therapy

Nasal high-flow (NHF) oxygen therapy is a technique that is increasingly used to treat ARF [12]. It has also recently been evaluated in the long-term domiciliary treatment of stable hypercapnic COPD patients. It was shown to improve health-related quality of life and to lower PCO₂ values to a certain extent (adjusted treatment effect, -4.1 mm Hg, 95% confidence interval, -6.5 to -1.7 mm Hg) [13]. tcPCO₂ measurements were used to demonstrate a decrease in PCO₂ values in different NHF usage settings: short-term use in chronic stable COPD patients [14], in tracheostomized COPD patients [15] and in stable COPD patients at night [16]. Therefore, tcPCO₂ measurements can be used to measure the extent of PCO₂ decline during NHF therapy.

Weaning

In order to wean a patient from invasive ventilation, spontaneous breathing trials (SBT) are performed on a regular basis [17]. In this context, it is crucial that PCO₂ is also monitored in addition to other parameters so that a decision can be made on when an SBT needs to be terminated. tcPCO₂ was evaluated in patients undergoing prolonged weaning in a trial by Schwarz et al [18]. It proved to be a suitable method for monitoring PCO₂ in patients undergoing invasive mechanical ventilation (MV) and prolonged weaning within the predefined clinically acceptable range of +/- 4 mmHg. This was in contrast to etCO₂, which significantly underestimated PaCO₂ (figure 2), particularly in patients with COPD due to a ventilation perfusion mismatch (mean difference -9.0 mmHg, limits of agreement -17.2 to -0.8 mmHg) [18].




Figure 2: Patient with COPD, PaCO₂ (black boxes), tcPCO₂ (orange line) and etCO₂ (dashed grey line) [18]

Outpatient setting and home care environment

Internationally, there is a trend towards a greater focus on outpatient care of patients with home mechanical ventilation (HMV). This relates to setting up MV at home rather than in the hospital [19] as well as outpatient check-ups of HMV [20]. In the Netherlands, tcPCO₂ has become the standard method for monitoring ventilation in the home environment [19] and is therefore able to replace blood gas analysis for stable patients receiving noninvasive HMV. In Germany, check-ups of noninvasive HMV in an inpatient setting still represent the standard of care [21], but there are initial pilot projects in place aimed at following up patients in an outpatient center. In this context, tcPCO₂ is used to evaluate short-term changes following adjustment of the ventilation settings [20]. As public healthcare systems are transitioning towards outpatient rather than inpatient care, noninvasive, easy-to-use, telemonitoring-compatible devices for monitoring gas exchange will be very useful in the future.

Functional assessment

tcPCO₂ was also shown to reliably monitor PCO₂ continuously in very severe COPD patients during a 6-minute walk test (6MWT); this could be of predictive value for the course of the disease and therefore provide useful information on the severity of respiratory pump insufficiency [22]. Interestingly, the trial revealed a very heterogeneous PCO₂ response during the 6MWT: While 24% of the patients preserved their PCO₂ values, 26% had reduced PCO₂ due to hyperventilation and 50% had increased PCO₂ [22].

Sedation during endoscopies

Numerous trials have demonstrated the use of tcPCO₂ for monitoring patients at risk of hypercapnia undergoing sedation during an endoscopy [1]. It was shown to detect hypoventilation and hypercapnia during bronchoscopy [23-25], as well as during colonoscopies [26] and thoracoscopies [27].

Summary

tcPCO₂ has proven to be beneficial in numerous clinical applications:

- it is **able to reduce the frequency of blood gas analysis** and to show a continuous trend of PCO₂ between blood gas analyses
- it provides **short and long-term feedback** of changes in alveolar ventilation in mechanically ventilated patients (after adjustment of (noninvasive) ventilation as well as during episodes of spontaneous breathing in weaning)
- it can be used to **measure the extent of PCO₂ decline during NHF therapy**
- it is able to **detect hypoventilation or hyperventilation** during clinical exercise tests (e.g. 6MWT)
- it is **superior to etCO₂ in patients with lung diseases and ventilation perfusion mismatch**
- it is reliable in monitoring PCO₂ in an **outpatient setting or in a home environment**, where blood gas analysis is not feasible

The SenTec Digital Monitoring System has proven to be accurate, reliable, easy to use and patient-friendly.

Literature

1. Huttmann SE, Windisch W, Storre JH: **Techniques for the measurement and monitoring of carbon dioxide in the blood.** *Ann Am Thorac Soc* 2014, 11(4): 645-652.
2. Nardi J, Prigent H, Adala A, Bohic M, Lebagry F, Quera-Salva MA, Orlikowski D, Lofaso F: **Nocturnal oximetry and transcutaneous carbon dioxide in home-ventilated neuromuscular patients.** *Respir Care* 2012, 57(9):1425-1430.
3. Kesten S, Chapman KR, Rebeck AS: **Response characteristics of a dual transcutaneous oxygen/carbon dioxide monitoring system.** *Chest* 1991, 99(5): 1211-1215.
4. Freeman ML, Hennessy JT, Cass OW, Pheley AM: **Carbon dioxide retention and oxygen desaturation during gastrointestinal endoscopy.** *Gastroenterology* 1993, 105(2):331-339.
5. Vargo JJ, Zuccaro G, Jr., Dumot JA, Conwell DL, Morrow JB, Shay SS: **Automated graphic assessment of respiratory activity is superior to pulse oximetry and visual assessment for the detection of early respiratory depression during therapeutic upper endoscopy.** *Gastrointestinal endoscopy* 2002, 55(7):826-831.
6. Windisch W, Waltersbacher S, Simon K, Geiseler J, Sitter H, German Society for P: **Guidelines for non-invasive and invasive mechanical ventilation for treatment of chronic respiratory failure. Published by the German Society for Pneumology (DGP).** *Pneumologie* 2010, 64(10):640-652.
7. Storre JH, Magnet FS, Dreher M, Windisch W: **Transcutaneous monitoring as a replacement for arterial PCO₂ monitoring during nocturnal non-invasive ventilation.** *Respir Med* 2011, 105(1):143-150.
8. Cox M, Kemp R, Anwar S, Athey V, Aung T, Moloney ED: **Non-invasive monitoring of CO₂ levels in patients using NIV for AECOPD.** *Thorax* 2006, 61(4):363-364.
9. Storre JH, Steurer B, Kabitz HJ, Dreher M, Windisch W: **Transcutaneous PCO₂ monitoring during initiation of noninvasive ventilation.** *Chest* 2007, 132(6): 1810-1816.
10. Lermuzeaux M, Meric H, Sauneuf B, Girard S, Normand H, Lofaso F, Terzi N: **Superiority of transcutaneous CO₂ over end-tidal CO₂ measurement for monitoring respiratory failure in nonintubated patients: A pilot study.** *Journal of critical care* 2016, 31(1):150-156.
11. Bobbia X, Claret PG, Palmier L, Robert M, Granpierre RG, Roger C, Yan J, Ray P, Seb-bane M, Muller L et al: **Erratum: Concordance and limits between transcutaneous and arterial carbon dioxide pressure in emergency department patients with acute respiratory failure: a single-center, prospective, and observational study.** *Scandinavian journal of trauma, resuscitation and emergency medicine* 2015, 23:77.
12. Zhu Y, Yin H, Zhang R, Wei J: **High-flow nasal cannula oxygen therapy versus conventional oxygen therapy in patients with acute respiratory failure: a systematic review and meta-analysis of randomized controlled trials.** *BMC Pulm Med* 2017, 17(1):201.
13. Nagata K, Kikuchi T, Horie T, Shiraki A, Kitajima T, Kadowaki T, Tokioka F, Chohnabayashi N, Watanabe A, Sato S et al: **Domiciliary High-Flow Nasal Cannula Oxygen Therapy for Patients with Stable Hypercapnic Chronic Obstructive Pulmonary Disease. A Multicenter Randomized Crossover Trial.** *Ann Am Thorac Soc* 2018, 15(4):432-439.
14. Fraser JF, Spooner AJ, Dunster KR, Anstey CM, Corley A: **Nasal high flow oxygen therapy in patients with COPD reduces respiratory rate and tissue carbon dioxide while increasing tidal and end-expiratory lung volumes: a randomised crossover trial.** *Thorax* 2016, 71(8):759-761.

15. Fricke K, Tatkov S, Domanski U, Franke KJ, Nilius G, Schneider H: **Nasal high flow reduces hypercapnia by clearance of anatomical dead space in a COPD patient.** *Respiratory medicine case reports* 2016, 19:115-117.
16. Biselli PJ, Kirkness JP, Grote L, Fricke K, Schwartz AR, Smith P, Schneider H: **Nasal high-flow therapy reduces work of breathing compared with oxygen during sleep in COPD and smoking controls: a prospective observational study.** *J Appl Physiol (1985)* 2017, 122(1):82-88.
17. Boles JM, Bion J, Connors A, Herridge M, Marsh B, Melot C, Pearl R, Silverman H, Stanchina M, Vieillard-Baron A et al: **Weaning from mechanical ventilation.** *Eur Respir J* 2007, 29(5): 1033-1056.
18. Schwarz SB, Windisch W, Magnet FS, Schmoor C, Karagiannidis C, Callegari J, Huttmann SE, Storre JH: **Continuous non-invasive PCO₂ monitoring in weaning patients: Transcutaneous is advantageous over end-tidal PCO₂.** *Respirology (Carlton, Vic)* 2017, 22(8):1579-1584.
19. Hazenberg A, Kerstjens HA, Prins SC, Vermeulen KM, Wijkstra PJ: **Initiation of home mechanical ventilation at home: a randomised controlled trial of efficacy, feasibility and costs.** *Respir Med* 2014, 108(9):1387-1395.
20. Schwarz SB, Callegari J, Hamm C, Windisch W, Magnet FS: **Is Outpatient Control of Long-Term Non-Invasive Ventilation Feasible in Chronic Obstructive Pulmonary Disease Patients?** *Respiration* 2017.
21. Windisch W, Dreher M, Geiseler J, Siemon K, Brambring J, Dellweg D, Grolle B, Hirschfeld S, Kohnlein T, Mellies U et al: **[Guidelines for Non-Invasive and Invasive Home Mechanical Ventilation for Treatment of Chronic Respiratory Failure - Update 2017].** *Pneumologie* 2017, 71(11):722-795.
22. Andrianopoulos V, Vanfleteren LE, Jarosch I, Gloeckl R, Schneeberger T, Wouters EF, Spruit MA, Kenn K: **Transcutaneous carbon-dioxide partial pressure trends during six-minute walk test in patients with very severe COPD.** *Respir Physiol Neurobiol* 2016, 233:52-59.
23. Chhajed PN, Rajasekaran R, Kaegi B, Chhajed TP, Pflimlin E, Leuppi J, Tamm M: **Measurement of combined oximetry and cutaneous capnography during flexible bronchoscopy.** *Eur Respir J* 2006, 28(2): 386-390.
24. Evans EN, Ganeshalingam K, Ebden P: **Changes in oxygen saturation and transcutaneous carbon dioxide and oxygen levels in patients undergoing fiberoptic bronchoscopy.** *Respir Med* 1998, 92(5):739-742.
25. Dreher M, Ekkernkamp E, Storre JH, Kabitz HJ, Windisch W: **Sedation during flexible bronchoscopy in patients with pre-existing respiratory failure: Midazolam versus Midazolam plus Alfentanil.** *Respiration* 2010, 79(4):307-314.
26. Heuss LT, Chhajed PN, Schnieper P, Hirt T, Beglinger C: **Combined pulse oximetry/cutaneous carbon dioxide tension monitoring during colonoscopies: pilot study with a smart ear clip.** *Digestion* 2004, 70(3):152-158.
27. Chhajed PN, Kaegi B, Rajasekaran R, Tamm M: **Detection of hypoventilation during thoracoscopy: combined cutaneous carbon dioxide tension and oximetry monitoring with a new digital sensor.** *Chest* 2005, 127(2):585-588.

sentec

SenTec AG

Ringstrasse 39

4106 Therwil

Switzerland

www.sentec.com

 Made in Switzerland

