

Vital patient data at a glance



Extend your CIDELEC PSG System by Sentec – **tcPCO_2 , drift corrected tcPCO_2 , SpO_2 and Pulse Rate (PR)** – measurements easily synchronized with CIDELEC P(S)G readings and directly integrated into the CIDELEC software.






Sentec Connectivity CIDELEC Analysis Software

- Sentec Digital Monitor (SDM) readings are duplicated to CIDELEC Analysis Software
- Easy setup for perfect synchronization of the SDM and P(S)G readings by one software (CIDELEC)
- Good resolution of tcPCO_2 / SpO_2 signals
- Retrospective correction of residual tcPCO_2 drift implemented into the PSG Analysis Software

Care with
Confidence

Sentec Connectivity

to CIDELEC Analysis Software

Required Components	System Requirements	Description
Sentec Digital Monitor (SDM) 	Software version SMB V07.xx to V08.04.xx	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
CIDELEC Analysis Software 	Version 2.2.4 or higher (compatible with CID102L, CID102L8, CID102L8D, CID-LXe, CID-LXe-206d, CID-LX, CID-LX-206d)	Synchronizes and integrates SDM readings of $tcPCO_2$, $tcPCO_2$ drift corrected, SpO_2 and Pulse Rate (PR) with CIDELEC P(S)G readings directly into the CIDELEC Software
Standard RS232 cable or USB2Serial converter cable 		Connection between SDM and CIDELEC Analysis Software

Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to connectivity@sentec.com or consult your CIDELEC representative.

Please visit www.sentec.com for more information.

Your local distributor:

sentec.

Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
 Made in Switzerland

Sentec Connectivity Support Team:
Tel +41 61 726 97 60
connectivity@sentec.com



20, rue des Métiers
49130 Sainte-Gemmes-sur-Loire
France
service.commercial@cidelec.net
www.cidelec.net
Tel +33 2 41 66 20 88

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.