


More vital patient data at a glance



SenTec Connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO₂**, **tcPO₂**, **SpO₂**, **Pulse Rate (PR)**, **Pulsation Index (PI)**, **Heating Power (HP)** - to the Dräger Integrated Care Manager (ICM).






Dräger

SenTec Connectivity

Dräger Integrated Care Manager (ICM)

- Displays and stores near real-time data from the SenTec Digital Monitor
- The SenTec Digital Monitor is automatically identified when connected to the Integrated Care Manager
- Integrated Care Manager supports and provides a dedicated and ergonomic clinical workflow for clinicians

Dräger Integrated Care Manager (ICM)

Required Components	System Requirements	Description
<p>SenTec Digital Monitor (SDM)</p> 	<p>Software version SMB V08.00/V08.01/V08.02</p>	<p>Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation</p>
<p>Integrated Care Manager (ICM)</p> 	<p>Version 9.2</p>	<p>ICM is a medical device software with Device Class IIa. For further detailed information contact your Dräger representative</p> <p>Displays and stores customizable data from SDM</p>
<p>Standard RS232 cable or USB2Serial converter cable</p> 		<p>Connects the SDM with the PC on which ICM is running</p>

SenTec Connectivity Support Team

We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to connectivity@sentec.com or consult with your Dräger representative.

Please visit www.sentec.com for more information.

sentec

Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
+ Made in Switzerland

Connectivity Support Team
+49 381 3677 96127
connectivity@sentec.com

Dräger

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck
Germany
www.draeger.com

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <https://www.sentec.com/info/legal-disclaimer/>