

## Vital patient data at a glance



Sentec connectivity enables remote interpretation of **tcPCO<sub>2</sub>** along with Trilogy Evo ventilation parameters on the Philips Care Orchestrator platform.



## Sentec Connectivity Philips Trilogy Evo

- Philips Trilogy Evo logs Sentec transcutaneous CO<sub>2</sub> monitoring data in standby mode or during ventilation
- Transfer of Sentec monitoring data to the Philips Care Orchestrator therapy management platform
- Sentec tcCO<sub>2</sub> data and Trilogy Evo ventilation data in the same report

# Sentec Connectivity

Philips Trilogy Evo

Required Components	System Requirements	Description
<b>Sentec Digital Monitor (SDM)</b> 	Software version SMB V08.00 and higher	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation  <b>Note: When using the Sentec Digital Monitor in a home care environment, the accessory Sentec Isolation Transformer is required.</b>
<b>Philips Trilogy Evo, Trilogy Evo O2 and Trilogy EV300</b> 	Software version 1.06.02.00 and higher	Portable noninvasive and invasive ventilator that provides support for a wide range of adult and pediatric patients
<b>Philips Isolated USB-A to DB9 cable (REF 1127405)</b> 		Both Trilogy Evo USB ports can be used to connect the SDM
<b>Male-to-male DB9 Serial Null Modem Cable (off the shelf product)</b>		

## Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult your Philips Respironics representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.

Your local distributor:

**sentec.**

Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)

 Made in Switzerland

Sentec Connectivity Support Team:  
Tel +41 61 726 97 85  
[connectivity@sentec.com](mailto:connectivity@sentec.com)

Philips Respironics Customer Service:  
Tel +1-724-387-4000  
Tel +1-800-345-6443 (toll-free)  
Tel +49 8152 93060 (international)

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.