

More vital patient data at a glance



SenTec Connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO<sub>2</sub>, SpO<sub>2</sub>, Pulse Rate (PR)** - to Vendom 30/40.



## SenTec Connectivity

Air Liquide Healthcare Vendom 30/40

- Automatically detects the SenTec Digital Monitor
- Displays and stores real-time data from SenTec Digital Monitor

**Air Liquide Healthcare Vendom 30/40**

Required Components	System Requirements	Description
<p><b>SenTec Digital Monitor (SDM)</b></p> 	<p><b>Software version SMB V07.01 to V08.02</b></p>	<p><b>Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation</b></p>
<p><b>Vendom 30/40</b></p> 	<p><b>Software version V2.05 or higher</b></p>	<p><b>Stores and displays tcPCO2, SpO2 and Pulse Rate (PR) from SDM</b></p>
<p><b>Serial/USB port adapter cable (ref. YR115800)</b></p> 		<p><b>Connection between Vendom 30/40 and SDM</b></p>

SenTec Connectivity Support Team


We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult with your Air Liquide Healthcare representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.



Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)  
 Made in Switzerland

Connectivity Support Team  
+49 381 3677 96127  
[connectivity@sentec.com](mailto:connectivity@sentec.com)



Parc de Haute Technologie 6  
Rue Georges Besse  
92182 Antony Cedex  
France  
+33 1 40 96 66 0  
[www.airliquide.com/healthcare](http://www.airliquide.com/healthcare)

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <https://www.sentec.com/info/legal-disclaimer/>