

More vital patient data at a glance



SenTec Connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO<sub>2</sub>**, **tcPO<sub>2</sub>**, **SpO<sub>2</sub>**, **Pulse Rate (PR)**, **Pulsation Index (PI)**, **Heating Power (HP)** - to the Dräger Integrated Care Manager (ICM).






**Dräger**

## SenTec Connectivity

### Dräger Integrated Care Manager (ICM)

- Displays and stores near real-time data from the SenTec Digital Monitor
- The SenTec Digital Monitor is automatically identified when connected to the Integrated Care Manager
- Integrated Care Manager supports and provides a dedicated and ergonomic clinical workflow for clinicians

## Dräger Integrated Care Manager (ICM)

Required Components	System Requirements	Description
<p><b>SenTec Digital Monitor (SDM)</b></p> 	<p><b>Software version SMB V08.00/V08.01/V08.02</b></p>	<p><b>Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation</b></p>
<p><b>Integrated Care Manager (ICM)</b></p> 	<p><b>Version 9.2</b></p>	<p><b>ICM is a medical device software with Device Class IIa. For further detailed information contact your Dräger representative</b></p> <p><b>Displays and stores customizable data from SDM</b></p>
<p><b>Standard RS232 cable or USB2Serial converter cable</b></p> 		<p><b>Connects the SDM with the PC on which ICM is running</b></p>

### SenTec Connectivity Support Team

We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult with your Dräger representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.

**sentec**

Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)  
+ Made in Switzerland

Connectivity Support Team  
+49 381 3677 96127  
[connectivity@sentec.com](mailto:connectivity@sentec.com)

**Dräger**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck  
Germany  
[www.draeger.com](http://www.draeger.com)

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <https://www.sentec.com/info/legal-disclaimer/>