

MINDRAY BENEVIEW SERIES / SENTEC CONNECTIVITY




- OXYGENATION AND VENTILATION MONITORING (tcPCO₂, tcPO₂, SpO₂, PR)
- DIRECTLY INTERFACED INTO BENEVIEW SERIES

Mindray BeneView Series
requires SW version 5.25 or later

Benelink Module
requires SW version 1.6 or later

SenTec Digital Monitoring System (SDMS)
providing noninvasive and continuous readings of tcPCO₂, tcPO₂, SpO₂, Pulse Rate, Heating Power and Sensor Temperature

Software requirement
SW version SMB 08.00 or later



BENEFITS

- Ease of use thanks to the plug and play feature
- Possibility to duplicate related alarms and alarm settings on Mindray BeneView Series patient monitors
- Data management via Mindray PDS gateway or eGateway in HL7 protocol

PLEASE CONTACT US

sentec

Ringstr. 39
4106 Therwil | Switzerland
info@sentec.ch
www.sentec.ch
Tel.: +41-61-7269760
Fax: +41-61-7269761

mindray

healthcare within reach

Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Mindray Building, Keji 12th Road South, High-tech industrial park
Nanshan, Shenzhen 518057, P.R.China
service@mindray.com
www.mindray.com
Tel.: +86 755 81888998
Fax: +86 755 26582680

