

Vital patient data at a glance



Sentec connectivity enables the transmission of oxygenation and ventilation monitoring – **tcPCO₂, tcPO₂, SpO₂, Pulse Rate (PR)** – to Mindray BeneView (T5, T6, T8, T9) and BeneVision (N12, N15, N17, N19, N22) patient monitors.




Sentec Connectivity

Mindray BeneView & BeneVision series

- Sentec Digital Monitor readings are duplicated on BeneView & BeneVision patient monitors
- Possibility to duplicate related alarms and alarm settings on Mindray BeneView & BeneVision patient monitors
- Data management via Mindray PDS gateway or eGateway in HL7 protocol

Sentec Connectivity

Mindray BeneView & BeneVision series

Required Components	System Requirements	Description
Sentec Digital Monitor (SDM) 	Software version SMB V08.00 to V08.03	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
BeneView series (T5, T6, T8, T9)	Software version 5.25 or higher	Duplicates and stores tcPCO ₂ , tcPO ₂ , SpO ₂ as well as Pulse Rate (PR) from SDM
BeneVision series (N12, N15, N17, N19, N22)	Software version V01-06	
BeneLink Module 	Software Version 1.6 or higher (for BeneView series) / 1.7 or higher (for BeneVision series)	The BeneLink Module is inserted into Mindray patient monitor module rack
ID Adapter 		Connects RS232 cable with LAN cable
Mindray type C serial port adapting cable 		Connects SDM and ID Adapter
Standard LAN cable 		Connects BeneLink Module and ID Adapter

Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to connectivity@sentec.com or consult your Mindray representative.

Please visit www.sentec.com for more information.

Your local distributor:

sentec.

Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
 Made in Switzerland

Sentec Connectivity Support Team:
Tel +41 61 726 97 60
connectivity@sentec.com

mindray

Shenzen Mindray Bio-Medical Electronics Co., Ltd.
Mindray Building, Keji 12th Road South, High-tech industrial park, Nanshan, Shenzhen 518057, P.R.China
+86 755 81888998
service@mindray.com
www.mindray.com

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.