

## Vital patient data at a glance



Sentec Connectivity enables the transmission of oxygenation and ventilation monitoring  
– tcPCO<sub>2</sub>, SpO<sub>2</sub>,  
Pulse Rate (PR) –  
to Vendom 30/40.



## Sentec Connectivity

### Air Liquide Healthcare Vendom 30/40

- Automatically detects the Sentec Digital Monitor
- Displays and stores real-time data from Sentec Digital Monitor

# Sentec Connectivity

Air Liquide Healthcare Vendom 30/40

Required Components	System Requirements	Description
<b>Sentec Digital Monitor (SDM)</b> 	Software version V07.01 to V08.03	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
<b>Vendom 30/40</b> 	Software version V2.05 or higher	Stores and displays tcPCO <sub>2</sub> , SpO <sub>2</sub> and Pulse Rate (PR) from SDM
<b>USB2 Serial converter cable</b> 		Connection between Vendom 30/40 and SDM

## Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult your Air Liquide Healthcare representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.

Your local distributor:

**sentec.**

Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)  
 Made in Switzerland

Sentec Connectivity Support Team:  
Tel +41 61 726 97 60  
[connectivity@sentec.com](mailto:connectivity@sentec.com)

 **Air Liquide**  
HEALTHCARE

Parc de Haute Technologie 6  
Rue Georges Besse  
92182 Antony Cedex  
France  
+33 1 40 96 66 0  
[www.airliquide.com/healthcare](http://www.airliquide.com/healthcare)

**Note:** This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.