

SenTec Connectivity to Dräger Integrated Care Manager

More vital patient data at a glance

Dräger Integrated Care Manager (ICM)

- Displays and stores near real-time data from SDM
- The SenTec Digital Monitor is automatically identified when connected to Integrated Care Manager
- Integrated Care Manager supports and provides a dedicated and ergonomic clinical workflow for clinicians



SenTec connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO₂**, **tcPO₂**, **SpO₂**, **Pulse Rate (PR)**, **Pulsation Index (PI)**, **Heating Power (HP)** - to Dräger Integrated Care Manager (ICM).



to Dräger Integrated Care Manager (ICM)

Required Components	System Requirements	Description
<p>SenTec Digital Monitor (SDM)</p> 	<p>Software version SMB V08.00/V08.01/V08.02</p>	<p>Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation</p>
<p>Integrated Care Manager (ICM)</p> 	<p>Version 9.2</p>	<p>ICM is a Medical device Software with Device Class IIa. For further detailed information contact your Dräger representative.</p> <p>Displays and stores customizable data from SDM.</p>
<p>Standard RS232 cable or USB2Serial converter cable</p> 		<p>Connects SDM with the PC where ICM is running</p>

SenTec Connectivity Support Team

We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to connectivity@sentec.com or consult with your Dräger representative.

Please visit www.sentec.com for more detailed technical and service information.



Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
Tel +41 61 726 97 60
Fax +41 61 726 97 61

Connectivity Support Team
+49 381 3677 96127
connectivity@sentec.com



Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck | Germany
www.draeger.com

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <http://www.sentec.com/int/info/legal-disclaimer/>