

## Vital patient data at a glance



Sentec Connectivity enables the transmission of ventilation monitoring – **tcPCO<sub>2</sub>** – to Vivo 45 / 45LS and Nippy 4+.




## Sentec Connectivity Breas ventilator series

- Easy follow-up of treatment efficacy by combining readings from Vivo 45 , Vivo 45LS and Nippy 4+ ventilators and Sentec Digital Monitor
- True plug and play
- Display of trend data and related alarms on Breas ventilators
- Ventilators designed for pediatric and adult use in home or hospital environment

Care with  
Confidence

# Sentec Connectivity

Breas ventilator series

Required Components	System Requirements	Description
<b>Sentec Digital Monitor (SDM)</b> 	Software version SMB V07.01 to V08.03	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
<b>Breas ventilator series</b> Vivo 45 	FW 1.0.33 or higher	Stores and displays tcPCO <sub>2</sub> and related alarms from SDM
Vivo 45LS Vivo 45LS (US) 	FW 3.0.44 or higher FW 5.0.01 or higher	
Nippy 4+ (UK only) 	FW 4.0.51 or higher	
<b>Serial adapter cable (ref. no. 006179)</b> 		Connection between Breas ventilators and SDM

## Sentec Connectivity Support Team

Let's connect.

If you need support, please contact your local Sentec representative, send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult your Breas representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.

Your local distributor:

**sentec.**

Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)  
 Made in Switzerland

Sentec Connectivity Support Team:  
Tel +41 61 726 97 60  
[connectivity@sentec.com](mailto:connectivity@sentec.com)

**BREAS**

Företagsvägen  
SE-435 33 Mölnlycke  
Sweden  
Phone +46 31 86 88 00  
[www.breas.com](http://www.breas.com)

**Note:** This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.