

More vital patient data at a glance



SenTec connectivity enables the transmission of oxygenation and ventilation monitoring - **tcPCO<sub>2</sub>, tcPO<sub>2</sub>, SpO<sub>2</sub>, Pulse Rate (PR) and Absolute Heating Power (AHP)** - to Philips patient monitors.









## SenTec Connectivity

### Philips patient monitor series

- SenTec Digital Monitor readings are duplicated on Philips patient monitors
- Related alarms and alarm settings are displayed on Philips patient monitors
- Possibility to duplicate the plethysmographic waveform

**Philips patient monitor series**

Required Components	System Requirements	Description
<b>SenTec Digital Monitor (SDM)</b> 	Software version SMB V07.01 to V08.02	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
<b>Patient monitors with IntelliBridge EC10 module rack</b> 	Software version Rev. H.15 or higher	Displays and stores tcPCO <sub>2</sub> , tcPO <sub>2</sub> , SpO <sub>2</sub> , Pulse Rate (PR) and Absolute Heating Power (AHP) on IntelliVue™ series
	All software versions	MP40/50/60/70/80/90 MX600/700/800
		MX400/450/500/550
<b>IntelliBridge EC10 Module</b> 		The IntelliBridge EC10 Module is inserted into the Philips patient monitor module rack
<b>IntelliBridge EC5 ID Module*</b> 		One module connects the SDM with the Patch cable
<b>IntelliBridge Adapter*</b> 		One adapter connects the SDM with the Patch cable
<b>Patch cable</b> 		Connection between IntelliBridge EC5 Module / IntelliBridge Adapter and IntelliBridge EC 10 Module

\* depending on the SDM motherboard, an EC5 Module or an IntelliBridge Adapter is needed. For further information, please refer to the VueLink / IntelliBridge Installation Manual.

SenTec Connectivity Support Team


We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to [connectivity@sentec.com](mailto:connectivity@sentec.com) or consult your Philips representative.

Please visit [www.sentec.com](http://www.sentec.com) for more information.



Ringstrasse 39  
4106 Therwil  
Switzerland  
[www.sentec.com](http://www.sentec.com)  
 Made in Switzerland

Connectivity Support Team  
Tel +49 381 3677 96127  
[connectivity@sentec.com](mailto:connectivity@sentec.com)

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <https://www.sentec.com/info/legal-disclaimer/>