

SenTec Connectivity to Cidelec

More vital patient data at a glance

Cidelec Analysis Software




- SenTec Digital Monitor readings are duplicated to Cidelec Analysis Software
- Easy setup for perfect synchronization of the SDM and P(S)G readings by one software (Cidelec)
- Good resolution of tcPCO₂/SpO₂ signals
- Retrospective correction of residual tcPCO₂ drift implemented into the PSG Analysis Software



Extend your Cidelec PSG System by SenTec - **tcPCO₂**, **drift corrected tcPCO₂**, **SpO₂** and **Pulse Rate (PR)** - measurements easily synchronized with Cidelec P(S)G readings and directly integrated into the Cidelec software.



to Cidelec Analysis Software

Required Components	System Requirements	Description
SenTec Digital Monitor (SDM) 	Software version SMB V07.xx to V08.01	Provides noninvasive and continuous real-time monitoring of ventilation and oxygenation
Cidelec Analysis Software 	Version 2.2.4 or higher (compatible with CID102L, CID102L8, CID102L8D, CID-LXe, CID-LXe-206d, CID-LX, CID-LX-206d)	Synchronizes and integrates SDM readings of tcPCO2, tcPCO2 drift corrected, SpO2 and Pulse Rate (PR) with Cidelec P(S)G readings directly into the Cidelec Software
Standard RS232 cable or USB2Serial converter cable 		Connection between SDM and Cidelec Analysis Software

SenTec Connectivity Support Team

We help you get connected

Your local distributor

If you need any support, please contact your local SenTec representative or send an e-mail to connectivity@sentec.com.

Please visit www.sentec.com for more detailed technical and service information.

sentec

Ringstrasse 39
4106 Therwil
Switzerland
www.sentec.com
Tel +41 61 726 97 60
Fax +41 61 726 97 61

Connectivity Support Team
Tel +49 381 3677 96127
connectivity@sentec.com



20, rue des Métiers
49130 Sainte-Gemmes-sur-Loire
France
service.commercial@cidelec.net
www.cidelec.net
Tel +33 2 41 66 20 88

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice. <http://www.sentec.com/int/info/legal-disclaimer/>