

SenTec warranty and license periods

REF	Product name	License period	Warranty period		Expected useful life ²
			After purchase by retailer	After delivery to end-customer ¹	
SDM	SenTec Digital Monitor	N.A.	24 months	24 months extension to 48 months possible	7 years
VS-A/P/N	V-Sign™ Sensor 2		21 months excluding the cable	18 months excluding the cable	up to 36 months
OV-A/P/N	OxiVenT™ Sensor	Life Time: 12 months ³ Usage Time: 4320 hours ³	15 months	12 months 4320 hours ³ excluding the cable	12 months
	Other accessories excluding cables	N.A.	24 months	24 months	7 years

¹ Date of invoice, installation or acceptance into service, whichever dates first.

² The average amount of time that an item is estimated to function when installed new, assuming correct/diligent handling and routine maintenance is practiced. For the avoidance of doubt, the indicated expected useful life time is for information purposes only and as such does not constitute, imply or establish any warranty or guarantee.

³ after first ex-factory use. The use of OxiVenT™ Sensor is granted in the sense of a license to the customer and its use is only granted and permitted during the Life Time and Usage Time (license period). With expiration of the Life Time the OxiVenT™ Sensor's operation automatically terminates. With expiration of the Usage Time the OxiVenT Sensor's PO2

measurement is automatically disabled. For the avoidance of doubt, the license is granted for the Life Time only and is not renewed or prolonged in case the OxiVenT™ Sensor is repaired or replaced by SenTec due to defect.

SenTec warranty and license periods

REF	Product name	Warranty period	Minimal shelf life ex. SenTec AG	Reusable
EC-MI	Ear Clip	Quality of material and assembly is warranted within the expiration date printed on the label	12 months	No
MAR-MI and MARE-MI	Multi-Site Attachment Ring for mature and intact skin		12 months	No
MAR-SF and MARE-SF	Multi-Site Attachment Ring for sensitive and fragile skin		9 months	No
SA-MAR	Staysite™ Adhesive for Multi-Site Attachment Ring		9 months	No
MC-R	Membrane Changer Reloadable		12 months	Yes, max. 10 times reloadable with MC-I
MC-I	Membrane Changer Insert		12 months	No
GEL-04	Contact Gel		12 months	Yes
GEL-SD	Single Dose Contact Gel		12 months	No
GAS-0812	Service Gas		12 months	Yes, for about one month ⁴

⁴ Depending on use scenarios and sensor condition.